



## Freezer Racks for Mechanical Freezers



Mechanical Freezer Racks

### Mechanical Freezer Racks

- Vertical racks for chest-type freezers are available in 4, 7, 9 or 12-shelf sizes to accommodate 25, 81 or 100 ampule boxes.
- Racks feature a unique locking rod to hold boxes securely. These racks allow quick retrieval of samples and fast return of rack into the freezer.
- Heavy-duty construction and unique unibody design prevents warping and bending due to temperature stress.
- Designed for use with standard 2" high cryogenic freezing boxes. No special boxes are required.
- Offers a unique inventory system for your mechanical freezer when used in conjunction with Nalgene® plastic CryoBoxes.

- Inventory system
- Choose from a variety of sizes
- Heavy-duty design minimizes damage to freezer interiors



Ampule Storage Boxes

### Ampule Storage Boxes

- The unique numbered grid system and one-way lid provides easy ampule identification.
- Durable plastic construction keeps ampules secure in and out of liquid nitrogen.
- Choice of ampule sizes and capacities maximize storage.

## ORDERING INFORMATION

Model #	# of Shelves	Overall Dimensions Inches (cm)			US List Price
		W	H	D	
<b>Racks for 81- and 100-ampule Boxes*</b>					
HR509X8A	4	5.42 (13.8)	9.0 (22.9)	5.56 (14.1)	\$76.00
HR509X6A	7	5.42 (13.8)	15.55 (39.5)	5.56 (14.1)	\$92.00
HR509X10A	9	5.42 (13.8)	19.92 (50.6)	5.56 (14.1)	\$111.00
HR509X3A	12	5.42 (13.8)	26.5 (67.3)	5.56 (14.1)	\$126.00
<b>Racks for 25-ampule Boxes*</b>					
HR509X7A	4	3.18 (8.1)	9.0 (22.9)	3.29 (8.4)	\$72.00
HR509X5A	7	3.18 (8.1)	15.55 (39.5)	3.29 (8.4)	\$82.00
HR509X9A	9	3.18 (8.1)	19.92 (50.6)	3.29 (8.4)	\$93.00
HR509X4A	12	3.18 (8.1)	26.5 (67.3)	3.29 (8.4)	\$107.00
<b>CryoBoxes and Vials</b>					
CS509X4	Nalgene System 81-ampule box accommodates most internally threaded 1.0 and 1.5 ml vials (Case of 40)				\$306.00
CS509X24	Nalgene System 100-ampule box accommodates most internally threaded 1.0 and 1.5 ml vials (Case of 10)				\$87.00
CS509X5	2 ml Nalgene 81-ampule box (Case of 20)				\$154.00
CS509X3	2 ml Nalgene 25-ampule box (Case of 80)				\$394.00
AY509X32	Nalgene System 100 Cryogenic Vials-1 ml (Case of 500)				\$198.00
AY509X33	Nalgene System 100 Cryogenic Vials-1.5 ml (Case of 500)				\$198.00
AY509X34	Nalgene System 100 Cryogenic Vial Holder (Package of 1)				\$25.00
AY509X35	Nalgene System 100 Cryogenic Vial Holder (Package of 4)				\$95.00

\* When ordering racks it is important to measure the chamber height of your freezer and choose the rack height which allows for a minimum of 1" clearance between the top of the rack and the freezer lid.