



Transport and Storage System



Arctic Express Dual Shippers

Ship Dry and Store in Liquid Nitrogen in the Same Vessel!

- Arctic Express™ Dual Shipper Storage System transports non-hazardous samples by common carrier
- Holding time of up to 26 days when used as a dry shipper; up to 125 days when used as a laboratory cryovessel

☆ 36 Month Vacuum Warranty

Product Description

- Family of six cryovessels can be used as dry shippers or laboratory cryovessels.
- Color-coded can and cane system allows for easy retrieval of storage sample.
- Smaller models come complete with reusable "mushroom" shipping case to prevent damage.
- Lockable lid prevents unauthorized entry.

PRODUCT SPECIFICATIONS

	Dual 10	Dual 19	Dual 28	Dual 32	Dual 46
Model #	CK50920	CK50921	CK50922	CK50923	CK50925
Liquid Nitrogen Capacity (Liters)	10	18.5	28	32	46
Adsorbed Liquid Nitrogen Capacity (Liters)	3	3	8	7	8
Holding Time -- Laboratory (Days)	42	125	50	100	74
Holding Time -- Shipper (Days)	21	21	24	30	21
Sample Capacity (1.2 ml Ampules)	288	396	1440	1260	2700
Vessel					
Height: Inches (cm)	21.6 (54.9)	25.7 (65.3)	22 (55.9)	26.6 (67.6)	26.5 (67.3)
Diameter: Inches (cm)	10.2 (25.9)	14.5 (36.8)	18.2 (46.2)	18.2 (46.2)	20 (50.8)
Neck Opening: Inches (cm)	2.0 (5.08)	2.2 (5.6)	3.8 (9.7)	3.5 (8.9)	5.0 (12.7)
Weight					
Shipping: Lb. (kg)	31.8 (14.4)	45.3 (20.5)	48.3 (21.9)	52 (23.5)	55 (24.9)

ORDERING INFORMATION

	Dual 10	Dual 19	Dual 28	Dual 32	Dual 46
Model #	CK50920	CK50921	CK50922	CK50923	CK50925
US List Price	\$1,200.00	\$1,300.00	\$1,450.00	\$1,395.00	\$1,695.00