

Thermolyne
Cryogenics



IATA-Shipper Transport System



Arctic Express™ IATA- Approved Shippers

IATA Approved

- Approved for air shipment of potentially infectious material
- Secondary security compartment
- ☆ 36 Month Vacuum Warranty

Product Description

- Unique design approved for worldwide shipment of infectious materials by UN and IATA.
- Two-week holding time insures that precious samples arrive intact.
- Lockable cover with internal security compartment provides protection and sample isolation.

PRODUCT SPECIFICATIONS

Model #	CY50920
Liquid Nitrogen Capacity (Liters)	5.9
Holding Time (Days)	14
Ample Capacity (2 ml)	350
Vessel	
Height: Inches (cm)	24 (61.0)
Diameter: Inches (cm)	15 (38.1)
Neck Opening: Inches (cm)	8.5 (21.6)
Weight	
Shipping: Lb. (kg)	29 (13.2)

ORDERING INFORMATION

Model #	CY50920
US List Price	\$3,495.00