

Thermolyne  
Cryogenics



## Bio-Cane™ Can and Cane Systems



Bio-Cane Storage Systems

### It's Practical and Economical!

- 5 sizes to choose from
- Priced right for your budget
- ☆ 24 Month Warranty, 5 Year Vacuum Warranty

#### Product Description

- Available in five sizes to fit the requirements and space limitations of research or clinical laboratories.
- Super vacuum insulation, durable construction provide years of trouble-free service.
- Secure locking hasp prevents unauthorized entry.
- Uniform temperature environment for specimen storage.
- Every vessel includes six cans to accommodate ampules.

#### Additional Bio-Cane 20, 34, and 47 Features

- Polycarbonate lid withstands the years of rugged use.
- Color coded canister system for quick and easy inventory identification.

#### PRODUCT SPECIFICATIONS

Model #	Height Inches (cm)	Outside Diameter Inches (cm)	Neck Tube Opening Inches (cm)	Liquid Nitrogen Capacity (Liters)	System Capacity Static Evaporation Rate (Liters/Day)	Static Holding Time (Days)	Canister Diameter Inches (cm)	Vials* (5 per cane/6 per cane)	
								2.0 ml	1.2 ml
CK509X1	21.6 (54.8)	10.2 (25.9)	2 (5.0)	11	0.15	73	1.5 (3.8)	150	180
CK509X2	25.7 (65.3)	14.5 (36.8)	2 (5.0)	20.5	0.09	227	1.5 (3.8)	150	180
CK509X3	26.6 (67.6)	18.2 (46.2)	3.5 (8.8)	34.8	0.18	193	2.8 (7.1)	600	720
CK509X4	26.5 (67.3)	20 (50.8)	5 (12.7)	47.4	0.39	121	4 (10.2)	1260	1512
CK509X5	27.3 (69.3)	22 (55.9)	6.25 (15.8)	70.7	0.60	117	5 (12.7)	1560	1872

\*Quantities are based on canes without sleeves.

#### ORDERING INFORMATION

Model #	Description	Shipping Weight Lb. (kg)	Accessory Cart Model #	US List Price
CK509X1	Bio-Cane 11	16 (7.3)	N/A	\$897.00
CK509X2	Bio-Cane 20	29 (13.2)	AY509X9	\$917.00
CK509X3	Bio-Cane 34	41 (18.6)	AY509X9	\$1,353.00
CK509X4	Bio-Cane 47	47 (21.3)	AY509X9	\$1,473.00
CK509X5	Bio-cane 70	72 (32.7)	AY509X1	\$2,086.00