



## Bio-Cane™ Accessories



CryoClaw

### CryoCane™ (A)

- Aluminum canes hold five 2.0 ml or six 1.2 ml standard cryovials for storage in Bio-Cane systems.

### Cryogenic Vial Color Coders (B)

- Colorful polystyrene disks fit into top of standard Nalgene cryovials (1.2 or 2.0 ml).
- Provides quick identification of vials. Flat surface permits writing on top of vial.

### CryoCane Coder (C)

- Unique, colorful aluminum tabs keep track of CryoCane inventory.
- Flat writing surface for additional identification.
- Conveniently packaged in quantities of 100.

### CryoSleeve™ (D)

- Reusable polyvinyl chloride sleeve encloses a CryoCane for extra security.
- Will not become brittle during freezing or thawing.

### Cryogenic Vials (E)

- Polypropylene vials have high density polyethylene closures which outperform other vials for storage in vapor phase liquid nitrogen.

- Screw closure has internal thread; vial has external thread, reducing potential contamination.
- Deep-skirted closure allows single-handed aseptic technique while minimizing exposure to vial contents.
- White marking area and fill line.
- Vials and closures preassembled and radiation-sterilized.
- Non-cytotoxic and non-pyrogenic.

### CryoClaw™

- Removes ampules and boxes that have fallen from racks and cans in cryogenic vessels.
- Removes ampules from vessels having necks larger than two inches in diameter.
- Used for numerous applications around the lab from reaching overhead items to retrieving fallen objects between benches and behind instruments.
- Safely move radio-isotopes without touching the container.

## ORDERING INFORMATION

Model #	Description	US List Price	Model #	Description	US List Price
AY509X20	CryoCane(A), case/48, 1.2 and 2.0 ml capacity, 286 mm length	\$48.00	AY509X24	CryoCane Coders(C), case/100, Blue	\$25.00
AY509X27	Cryogenic Vials Closure Color Coders (B), pkg/100, White	\$14.00	AY509X25	CryoCane Coders(C), case/100, Green	\$25.00
AY509X28	Cryogenic Vials Closure Color Coders (B), pkg/100, Yellow	\$14.00	AY509X26	CryoCane Coders(C), case/100, Red	\$25.00
AY509X29	Cryogenic Vials Closure Color Coders (B), pkg/100, Blue	\$14.00	AY509X21	CryoSleeve(D), Case/100, 273 mm length	\$67.00
AY509X30	Cryogenic Vials Closure Color Coders (B), pkg/100, Green	\$14.00	AY509X13	Cryogenic Vials(E), case/500, 1.2 ml, 13.5 mm O.D., 38.1mm H	\$201.00
AY509X31	Cryogenic Vials Closure Color Coders (B), pkg/100, Red	\$14.00	AY509X12	Cryogenic Vials(E), case/500, 2.0 ml, 13.5 mm O.D., 48.3 mm H	\$201.00
AY509X22	CryoCane Coders(C), case/100, White	\$25.00	AY509X18	CryoClaw	\$115.00
AY509X23	CryoCane Coders(C), case/100, Yellow	\$25.00	AY509X19	Replacement suction cups for CryoClaw	\$14.00