



Locator® Rack and Box Systems



Keep Track of Your Samples—Easily!

- Indexed rack and box system
- Stores up to 3600 samples
- ☆ 36 Month Warranty

Locator Storage Systems

Product Description

- Locator systems available in three sizes with capacities up to 3600 vials.
- Unique hanging rack and gridded box design allows tracking and easy retrieval of samples.
- Suited for both manual and computerized inventory record-keeping methods.
- Advanced vacuum insulation minimizes liquid nitrogen evaporation and reduces operating costs.
- Unbeatable temperature uniformity with advanced vacuum insulation, assures all samples are stored below -190°C. (Even when less than 2" (5 cm) of liquid nitrogen remains in the vessel.)
- Backed by Thermolyne's exclusive three-year warranty, Locator's durable construction provides years of trouble-free service.
- Compact size easily fits into any lab.
- Secure locking hasp prevents unauthorized entry.
- An optional liquid nitrogen level monitor safeguards your valuable samples.

PRODUCT SPECIFICATIONS

	Locator Jr.	Locator 4	Locator 8		Locator Jr.	Locator 4	Locator 8
Height: Inches (cm)	28.9 (73.3)	37.5 (95.3)	38.3 (97.2)	Racks Included With Vessel (2-inch box rack)	4	4	8
Diameter: Inches (cm)	22 (55.8)	22 (55.8)	22 (55.8)	Boxes Per Rack	4	9	9
Neck Diameter: Inches (cm)	8.5 (21.5)	8.5 (21.5)	6 (15.2)	Box Capacity (Ampules)	100	100	25
Working Liquid Nitrogen Capacity (Liters)	51	111	111	Box Dimensions Inches (cm)	5.25 x 5.25 (13.3 x 13.3)	5.25 x 5.25 (13.3 x 13.3)	3 x 3 (7.6 x 7.6)
Average Static Liquid Nitrogen Consumption Rate (Liters/Day)	0.85	0.99	0.60	Boxes Needed For System	16	36	72
Static Holding Time (Days)	60	112	185	*Shipping Weight: Lb. (kg)	90 (40.8)	115 (52.2)	117 (53.1)
System Capacity (Ampules)	1600	3600	1800	*Add 4 lbs. for units with level monitor.			

ORDERING INFORMATION

Model #	Description	US List Price
CY50925	Locator Jr. with 2-inch box racks	\$2,709.00
CY50935	Locator 4 with 2-inch box racks	\$3,193.00
CY50945	Locator 8 with 2-inch box racks	\$3,058.00
CY50955	Locator Jr. with level monitor/120 Volt	\$3,729.00
CY50965	Locator 4 with level monitor/120 Volt	\$4,254.00
CY50975	Locator 8 with level monitor/120 Volt	\$4,095.00
AY509X1	Wheeled Accessory Cart (5 inches high) For all Locators (except Locator 6 Plus - AY509X1-70) and Biocane 70	\$216.00

240 Volt CE marked Level Monitors are available. Contact Barnstead International.