



## Locator® Accessories



### Level Monitor

- Safeguard irreplaceable biological specimens with the most advanced liquid nitrogen level monitor available.
- Unlike traditional monitoring devices, this ultrasonic device eliminates contact with liquid nitrogen reducing liquid nitrogen loss due to conduction.
- LED level indicator panel continuously displays liquid nitrogen volume in 1/8 increments.
- Remotable alarm notifies you of low liquid nitrogen levels even at times when no one is in the building.
- Audible and visual alarms provide early warning of low level conditions. Indicator light flashes when approximately 2" (5 cm) of liquid nitrogen remains in the vessel. An audible alarm sounds in unison with flashing low-level indicator light when level drops to about 1" (2.5 cm).



### Locator Log

- Designed for quick, reliable sample identification storage and relocation.
- Three section format for Locator and Locator Plus systems recordkeeping.
  - Section 1: Serves as a helpful reminder for refilling vessels.
  - Section 2: Locates your samples by identification keys.
  - Section 3: Inventories samples and pinpoints exact location.



### CryoClaw™

- Designed to remove ampules and boxes that have fallen from racks and cans in cryogenic vessels having necks larger than two inches in diameter.
- Reach overhead items or retrieve fallen objects between benches and behind instruments.



### ACCESSORIES

Model #	Description	US List Price	Model #	Description	US List Price
AY509X1	Wheeled Accessory Cart, 5 inches high (12.7 cm) For all Locators (except Locator 6 Plus - AY509X1-70) and Biocane 70	\$216.00	HR509X20A	Locator 4, 4-inch box rack	\$139.00
CS509X24	Nalgene® System™ 100-ampule box accommodates most internally threaded 1.0 and 1.5 ml vials (Case of 10) Size: 5.25 x 5.25 x 2.06 (13.34 X 13.34 X 5.23 cm)	\$87.00	LT509X5	Locator 4, Jr. & 4, 6 or Jr. Plus Logbook	\$149.00
CS509X3	2 ml Nalgene 25-ampule box (Case of 80) -- 2" box	\$394.00	LT509X6	Locator 8 & 8 Plus Logbook	\$161.00
CS509X4	2 ml Nalgene 81-ampule box (Case of 40) -- 2" box	\$306.00	AY509X18	CryoClaw	\$115.00
CS509X5	2 ml Nalgene 81-ampule box (Case of 20) -- 2" box	\$154.00	AY509X19	Replacement Suction Cups (2) for CryoClaw	\$14.00
CS509X10	5 ml Nalgene 81-ampule box (Each) -- 4" box	\$20.00	CN509X1	Level monitor for Locator Jr. (240 Volt)	\$1,061.00
AY509X32	Nalgene System 100 Cryogenic Vials--1 ml (Case of 500)	\$198.00	CN509X3	Level monitor for Locator Jr. (120 Volt)	\$1,061.00
AY509X33	Nalgene System 100 Cryogenic Vials--1.5 ml (Case of 500)	\$198.00	CN509X4	Level monitor for Locator 4 (240 Volt)	\$1,061.00
AY509X34	Nalgene SYstem 100 Cryogenic Vial Holder (Qty. 1)	\$25.00	CN509X6	Level monitor for Locator 4 (120 Volt)	\$1,076.00
AY509X35	Nalgene System 100 Cryogenic Vial Holder (Package of 4) Dimensions 4 x 7.75 x 0.8 (10.16 X 19.7 X 2.2 cm)	\$95.00	CN509X7	Level monitor for Locator 8 (240 Volt)	\$1,076.00
HR509X19A	Locator Jr. 4-inch box rack	\$110.00	CN509X9	Level monitor for Locator 8 (120 Volt)	\$1,076.00
<sup>1</sup> CSA			CN509X3-70 <sup>1</sup>	Level monitor for Locator Jr Plus (120 Volt)	\$1,108.00
<sup>2</sup> CE			CN509X6-70	Level monitor for Locator 4 Plus (120 Volt)	\$1,067.00
			CN509X10-70 <sup>1,2</sup>	Level monitor for Locator 6 Plus (120 Volt)	\$1,076.00
			CN509X9-70	Level monitor for Locator 8 Plus (120 Volt)	\$1,067.00