

Barnstead Deionization



E-Pure™



EASYpure 4-holder

Product Description

- Resistivities up to 18.2 megohm-cm and TOC levels below 10 ppb.
- Virgin polypropylene water pathways prevent recontamination.
- Recirculation pump continuously recirculates water throughout the entire system to maintain water purity and reduce the quantity of waste water required for rinse-up.
- Wall mounted; completely assembled for easy installation.

Low on Cost...High on Quality!

- 1/4 turn quick release canisters
- High flow rates
- Lowest operating cost of any system on the market
- Larger cartridges mean increased cartridge life



- 1/4 turn quick release canisters allow for easy cartridge replacement.
- 0.2 micron absolute hollow fiber filter removes bacteria and particulates that could alter your results.
- Digital resistivity meter is automatically temperature compensated to 25°C, preventing fluctuating resistivity readings due to temperature changes.
- Flow rate up to 2 liters per/min.
- High capacity cartridges provide low operating cost.

ORDERING INFORMATION

Model #	Electrical (50/60 Hz)	Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)	Max Inlet Pressure (psig)	Inlet Conn (Inches)	Inlet Temp Range	Max Flow Rates	US List Price
		W	H	D						
E-pure 3-Holder										
D4631	120 VAC	29 (73.7)	28 (71.1)	7.5 (19.1)	50 (22.7)	100	1/4 NPTF	40°-120°F	2.0 LPM	\$2,930.00
D4632-33	240 VAC	29 (73.7)	28 (71.1)	7.5 (19.1)	50 (22.7)	100	1/4 NPTF	4°-49°C	1.8 LPM	\$2,930.00
E-pure 4-Holder										
D4641	120 VAC	35 (88.9)	28 (71.1)	7.5 (19.1)	56 (25.4)	100	1/4 NPTF	40°-120°F	2.0 LPM	\$3,089.00
D4642-33	240 VAC	35 (88.9)	28 (71.1)	7.5 (19.1)	56 (25.4)	100	1/4 NPTF	4°-49°C	1.8 LPM	\$3,089.00

-33 Models CE Marked

ACCESSORIES

Cartridge Kits for E-pure Systems - Tap Feed (Final Filters Sold Separately)

Model #	Description	Contents	US List Price
D5029	3 Module Type I Cartridge Kit	1-D0835 Pretreatment Cartridge, 2-D5027 Ultrapure Cartridges	\$273.00
D5028	4 Module Type I Cartridge Kit	1-D0835 Pretreatment Cartridge, 1-D0803 High Capacity Cartridge, 2-D5027 Ultrapure Cartridges	\$354.00
D5022	3-Module Type I ORGANICfree Cartridge Kit	1-D0836 MACROpure Cartridge, 1-D5027 Ultrapure Cartridge, 1-D5021 ORGANICfree Cartridge	\$415.00
D5023	4 Module Type I ORGANICfree Cartridge Kit	1-D0836 MACROpure Cartridge, 1-D0803 High Capacity Cartridge, 1-D5027 Ultrapure Cartridge, 1-D5021 ORGANICfree Cartridge	\$418.00

Cartridge Kits for E-pure Systems - Pretreatment Feed (Final Filters sold separately)

D50227	4 Module Type I Cartridge Kit for Pretreated Feed	1-D0835 Pretreatment Cartridge, 2-D5027 and 1-D0809 Ultrapure Cartridge	\$360.00
D50228	4 Module Type I ORGANICfree Cartridge Kit for Pretreated Feed	1-D0836 MACROpure Cartridge, 1-D5027 Ultrapure Cartridges and 1-D0809 Ultrapure Cartridges, 1-D5021 ORGANICfree Cartridge	\$422.00
D3750	Final Filter		\$75.00