

Barnstead Deionization



Individual Cartridges

PRODUCT SPECIFICATIONS

Product Specifications for Pressurized Cartridges

Length Inches (cm)	Diameter Inches (cm)	Shipping Weight Lb. (kg)
17.2 (43.7)	3.4 (8.6)	5 (2.3)

Product Specifications for 1/2 Sized Cartridges

Length Inches (cm)	Diameter Inches (cm)	Shipping Weight Lb. (kg)
10.9 (27.7)	3.4 (8.6)	3 (1.4)

ORDERING INFORMATION

Model #	Description	Resin Type(s)	Where Used	US List Price
Pressurized Cartridges				
D0803	High capacity two-bed	Two bed, high capacity,	NANOpure®, E-pure™, B-pure™ Bantam, cooling loops	\$83.00
D0809	Ultrapure	Mixed bed	NANOpure II, older style NANOpure B-pure, Bantam, cooling loops	\$87.00
D0811	Oxygen removal	Oxygen	B-pure, cooling loops	\$90.00
D0813	Organic removal	Activated carbon	B-pure, cooling loops, Bantam	\$72.00
D0815	Cation removal	Cation	B-pure, cooling loops, Bantam	\$65.00
D0820	Organic free	Organic free	NANOpure II, older NANOpure systems	\$104.00
D5021	Organic free	Organic free	Quad-style NANOpure, E-pure	\$112.00
D0832	Pretreatment cartridge	Mixed bed, activated carbon	B-pure, still pretreatment	\$80.00
D0835	Combination colloid, organic	Macroreticular, activated carbon	NANOpure, E-pure, B-pure	\$82.00
D0836	Macropure	Macroreticular, activated carbon	NANOpure, E-pure, pretreatment	\$119.00
D8809	Two bed combination	Ultrapure, oxygen removal	B-pure, cooling loops	\$89.00
D8811	Three bed combination	Ultrapure, oxygen removal, organic removal	B-pure, cooling loops	\$88.00
D0760	Anion removal	Anion	B-pure, cooling loops, Bantam	\$156.00
D5027	Ultrapure	Mixed bed	Quad-style NANOpure, E-pure, ROpure ST	\$87.00
<i>*Cartridges feature all polypropylene construction and are designed to fit any NANOpure, E-pure, B-pure, PCS, Bantam holder or cooling loop system.</i>				
1/2 Size Cartridges for 1/2 Size B-Pure				
D50213	1/2 size Ultrapure	Mixed bed	1/2 B-pure	\$74.00
D50214	1/2 size oxygen removal	Oxygen removal	1/2 B-pure	\$84.00
D50215	1/2 size organic removal	Activated carbon	1/2 B-pure	\$61.00
D50216	1/2 size comb. oxygen and demineralization	Oxygen removal Mixed bed	1/2 B-pure	\$81.00
D50217	1/2 size comb. carbon and demineralization	Mixed bed Activated carbon	1/2 B-pure	\$79.00
NANOpure Infinity Individual Cartridges				
D50250	Pretreatment Cartridge R/O & Dist Feed	Activated Carbon	NANOpure Infinity	\$66.00
D50251	Pretreatment Cartridge DI & Tap Feed	Activated Carbon & Macroreticular Resin	NANOpure Infinity	\$85.00
D50252	Organic Free Cartridges	Organic Free Media	NANOpure Infinity	\$105.00
D50253	Ultrapure	Mixed Bed	NANOpure Infinity	\$112.00