



Replacement Cartridges



PUREpack Cartridge

Replacement Cartridges for Millipore® Water Systems

- Guaranteed performance
- Low price
- Available through your local laboratory distributor

PUREpack Deionization Replacement Cartridge

Replacements for Millipore water systems including: Milli-Q RG Water System (ZD51 115 84), Milli-Q Plus Water System (ZD52 115 84), Alpha-Q Water System (ZFHQ 060 01), Milli-Q UV Plus Water System (ZD60 115 UV), Milli-Q TOC Plus (SL1R 096 M4), Milli-Q PF Plus (ZD54 115 95), Milli-Q UF Plus (ZD53 115 95)

ORDERING INFORMATION

PUREpack Model #	Description	US List Price
D50265	PUREpack cartridge for reverse osmosis, distillation or ultrafiltration feedwater source (Includes final filter) Replacement for Millipore OPAK Cartridge (Millipore Catalog # CPMQ K05 R1.)	\$385.00
D50267	PUREpack cartridge for reverse osmosis, distillation or ultrafiltration feedwater source Replacement for Millipore OPAK Cartridge (Millipore Catalog # CPMQ 004 R1.)	\$350.00
D50266	PUREpack cartridge for deionized feedwater source (Includes final filter) Replacement for Millipore OPAK Cartridge (Millipore Catalog # CPMQ K05 D2.)	\$385.00
D50268	PUREpack cartridge for deionized feedwater source. Replacement for Millipore OPAK Cartridge (Millipore Catalog # CPMQ 004 D2.)	\$350.00

Barnstead Deionization



PUREpack RO Replacement Cartridge

Replacements for Millipore Reverse Osmosis systems including: Milli-RX® 20, Milli-RX 45, Milli-RX 75, Milli-RO 3 Plus, Milli-RO® 5 Plus, Milli-RO 8 Plus, Milli-RO 12 Plus, Milli-RO 30 Plus, Milli-RO 60 Plus, Milli-RO 90 Plus, AFS 10, AFS 18, AFS 30, AFS 200, AFS 300, AFS 1500, AFS 2000, AFS 3500, AFS 5500 RX

ORDERING INFORMATION

PUREpack Model #	Description	US List Price
D50269	PUREpack RO cartridge for compact systems and AFS analyzer feed systems Replacement for Millipore Catalog # CPRO P04 02.	\$282.00

PRODUCT SPECIFICATIONS

	Barnstead PUREpack1 w/filters	Barnstead PUREpack1	Barnstead PUREpack2 w/filters	Barnstead PUREpack2	Barnstead PUREpack RO
Cartridge Catalog Number	D50265	D50267	D50266	D50268	D50269
Final Filter Included	Yes	No	Yes	No	No
Grain Capacity/pak (CaCo3)	1121	1121	594	594	NA
Cost Per Liter(based on list price)	\$0.186	\$0.169	\$0.35	\$0.319	NA
Liters Produced (to 10 megohm endpoint)	2070	2070	1098	1098	NA
US List Price	\$385.00	\$350.00	\$385.00	\$350.00	\$282.00

ORDERING INFORMATION

PUREpack Model #	Description	US List Price
D50224	Replacement ultrapure cartridge kit compatible with Millipore 3 Bowl Milli-Q and 4 Bowl Milli-Q Plus Life Science systems. Kit includes two deionization cartridges, one carbon cartridge, and one final filter. Replacement for Millipore Catalog # CDMF 012 04 or CFUF 012 04.	\$193.00
D50225	Replacement ultrapure cartridge kit compatible with Millipore 4 Bowl Milli-Q for analytical and 5 Bowl Milli-Q Plus for life science/analytical.* Kit includes two deionization cartridges, one carbon cartridge, one prefilter and one final filter. Replacement for Millipore Catalog # CDMF 012 05 or CFIF 012 05.	\$266.00
D50226	Replacement ultrapure cartridge kit compatible with Millipore 4 Bowl Milli-Q for analytical and 5 Bowl Milli-Q Plus life science/analytical.* Kit includes two deionization cartridges, one carbon cartridge, one low TOC cartridge and one final filter. Replacement for Millipore Catalog # CDOF 012 05 or CFOF 012 05.	\$266.00
D8832	Individual deionization cartridges, compatible with above Millipore systems. Replacement for Millipore Catalog # CPMB 012 02.	\$56.00
D8834	Individual low TOC cartridge, compatible with above Millipore systems. Replacement for Millipore Catalog # CDEX 012 01.	\$78.00
D8831	Individual carbon cartridge, compatible with above Millipore systems. Replacement for Millipore Catalog # CDFC 012 04.	\$45.00
D8830	Individual prefilter, compatible with above Millipore systems. Replacement for Millipore Catalog # CWSC 01T P3.	\$94.00
D3750	Individual final filter, compatible with above Millipore systems. Replacement for Millipore Catalog # MPGL 04S K2 or CW03 012 C3.	\$75.00

*Does not include Pyrogard® ultrafilter, must be ordered separately from Millipore.

Millipore, Milli-Q, Milli-RX, Milli-RO, Millipak, Milligard, Ion-Ex, Pyrogard, QPAK and ROPAK are registered trademarks of Millipore Corporation and Millipore Investment Holdings Limited.