



## REPIPET® Dispensers



Standard and Low Silhouette  
 Dispensers



Universal REPIPET  
 Three cap and washer sets  
 sizes 38-430,  
 33-430,  
 28-430

### Our Best Dispensers!

- Standard dispensers are excellent for use with acids and solvents, except for HF
- Low silhouette dispensers are safe for dispensing and storing acids, except for HF
- 1% accuracy and 0.1% reproducibility at full scale guaranteed
- ☆ Glassware carries a 50% lifetime guarantee against breakage

### Product Description

- Dispenser components are fully autoclavable.
- Full range adjustability, 100-Division Scale.
- Includes a magnifying indicator for easy-reading accuracy.
- Pyrex® and Teflon® fluorocarbon construction for safe use with all liquids except hydrofluoric solutions.
- Models ending in A and AU include both magnifying indicator and metal indicator.

### Universal REPIPET Dispensers

- Dispense from any screw-cap container with an inner diameter opening greater than 19 mm or any container with a 24/40 or 29/42 standard taper by adding an adapter for a ground-glass joint.
- Supplied with 3-cap and washer sets to fit your container: 28-430, 33-430 & 38-430.
- 45 mm cap & washer set optional. #CV830X1B
- 24" (61 cm) of Teflon tubing to fit your container.

### Low Silhouette REPIPET Dispensers

- Ideal for shelf storage, fits easily in the refrigerator.
- Square bottle reduces the risk of damage from tipping.
- A second opening in the Halar®-fluorocarbon cap allows you to add reagent without disturbing the dispenser.
- Supplied with 1000 ml square amberized bottle.

### Standard REPIPET Dispensers

- Perfect all purpose choice for dispensing all reagents except hydrofluoric acids, including concentrated acids, concentrated alkalis, chlorinated hydrocarbons, and organic solvents.
- Supplied with a 950 ml round amber bottle.

### ORDERING INFORMATION

Model #	Size (ml)	Bottle Size (ml)	US List Price	Model #	Size (ml)	Bottle Size (ml)	US List Price	Model #	Size (ml)	Bottle Size (ml)	US List Price								
<b>Universal REPIPET Dispensers</b>																			
3001/2AU	0.5	Customer Supplied	\$306.00	<b>Low Silhouette REPIPET Dispensers</b>															
3001AU	1.0	Customer Supplied	\$306.00	3001/2AL	0.5	1000	\$315.00	<b>Standard REPIPET Dispensers</b>											
3005AU	5.0	Customer Supplied	\$314.00	3001AL	1.0	1000	\$302.00	3001/2A	0.5	950	\$302.00								
3010AU	10.0	Customer Supplied	\$314.00	3005AL	5.0	1000	\$320.00	3001A	1.0	950	\$302.00								
3020AU	20.0	Customer Supplied	\$336.00	3010AL	10.0	1000	\$320.00	3005A	5.0	950	\$306.00								
3050AU	50.0	Customer Supplied	\$409.00	3020AL	20.0	1000	\$327.00	3010A	10.0	950	\$306.00								
3100AU	100.0	Customer Supplied	\$511.00	3050AL	50.0	1000	\$415.00	3020A	20.0	950	\$314.00								
								3050A	50.0	950	\$389.00								