



## REPIPET® Jr. Fixed and Adjustable Volume



Model 770200 Fixed Volume Dispenser



Model 7010 Adjustable Volume Dispenser

### Our Lowest Priced Dispensers!

- Fixed volume available in 13 sizes
- Adjustable models are economical, plastic and unbreakable
- Excellent choice for aqueous solutions and HF

**Includes FREE support stand**

#### Fixed Volume Dispensers

- Specially designed recirculating system for priming without loss of reagent.
- Choose any of 13 sizes.
- Guaranteed 2% accuracy and 5% reproducibility.
- Each model includes one amber glass bottle and one plastic container.
- Adaptable for use with any reagent bottle having a 20-400, 22-400, 24-400, 28-400 or 38-400 screw cap.
- 5 interchangeable cap/adapters sets for container flexibility.
- Twelve inches (30.5 cm) of Teflon intake tubing.

#### Adjustable Volume Dispensers

- Economical and durable high density polyethylene and polypropylene dispenser.
- Time saving, semi-automatic reloading feature for repetitive lab tests.
- Most economical dispenser on the market today.
- Not suitable for use with oils, hydrocarbons, iodine, petroleum-based chemicals, oxidizing acids, alicyclics, solvents or scintillation cocktails.

### ORDERING INFORMATION

Fixed Volume				Adjustable Volume			
Model #	Size (ml)	Container Size	US List Price	Model #	Size (ml)	Container Size	US List Price
770100	0.1 ml	4 oz. (120 ml)	\$98.00	70021	2.0 ml	32 oz. (960 ml)	\$97.00
770200	0.2 ml	4 oz. (120 ml)	\$98.00	7005	5.0 ml	32 oz. (960 ml)	\$97.00
770500	0.5 ml	4 oz. (120 ml)	\$98.00	7010	10.0 ml	32 oz. (960 ml)	\$97.00
770700	0.7 ml	4 oz. (120 ml)	\$98.00				
770900	0.9 ml	8 oz. (240 ml)	\$98.00				
7701	1.0 ml	8 oz. (240 ml)	\$98.00				
77011/2	1.5 ml	8 oz. (240 ml)	\$98.00				
7702	2.0 ml	8 oz. (240 ml)	\$98.00				
77021/2	2.5 ml	8 oz. (240 ml)	\$98.00				
7703	3.0 ml	8 oz. (240 ml)	\$98.00				
7704	4.0 ml	16 oz. (480 ml)	\$98.00				
7705	5.0 ml	16 oz. (480 ml)	\$98.00				
7710	10.0 ml	16 oz. (480 ml)	\$98.00				