



## REPIPET® Single Stroke Dilutors



Model S81001U Dilutor

### Compare the Price to Automatic Dilutors!

- A specialized instrument for diluting samples
- Dilutes with a single action
- Universal style adapts to your container
- \* Glassware carries a 50% lifetime guarantee against breakage

### Product Description

- As fast and precise as expensive automatic dilutors.
- Dilutes with a single action; aspirating and diluting with just one cycle.
- Coupled together by a magnifying indicator coupling; diluent plunger and aspirator plunger raise simultaneously.
- 100-Division Scale.
- Both models feature magnifying indicators for easy reading accuracy.
- Pyrex® glass and Teflon® construction make REPIPET dilutors safe for use with all reagents except hydrofluoric solutions.
- Low Silhouette model includes a low-profile, 1000 ml square amberized bottle with two-opening lid.
- Low silhouette models are not suitable for strong organic solvents as the two-opening, Halar® lids may distort.
- Universal REPIPET dilutor supplied with three cap and washer sets and 24 inches of Teflon tubing: 28-430, 33-430 and 38-430 for use on your container. Bottle not included.

### ORDERING INFORMATION

Model #	Dilutor	Aspirator	Aspirated Volumes	Bottle Size (ml)	US List Price
<b>Low Silhouette</b>					
S80101/2L	1.0 ml	0.5 ml	1:1	1000	\$589.00
S81001/2L	10.0 ml	0.75 ml	13:1	1000	\$589.00
S81001L	10.0 ml	4.3 ml	2.3:1	1000	\$589.00
<b>Universal</b>					
S80101/2U	1.0 ml	0.5 ml	1:1	Customer Supplied	\$568.00
S81001/2U	10.0 ml	0.75 ml	13:1	Customer Supplied	\$568.00
S81001U	10.0 ml	4.3 ml	2.3:1	Customer Supplied	\$568.00