



## Air Displacement Pipets



Model T705B Pipet

**Includes three  
 disposable tips and  
 three extension tips**

### Chemically Resistant!

- All plastic construction
- Spring return, syringe type operation
- Resistant to a wide variety of chemicals

### Product Description

- Sturdy all-plastic construction makes it tough.
- Calibrated in 0.1 ml divisions, this macro transfer pipet is adjustable the full length of the scale.
- Accuracy within  $\pm 2\%$  at full volume. Reproducibility  $\pm 0.5\%$  at full volume.
- Not to be used with oils, hydrocarbons, toluene, petroleum based chemicals, oxidizing and concentrated acids, aromatic solvents, alicyclics, mercaptol reagents or iodine.
- Supplied with three disposable tips and three extension tips.

### ORDERING INFORMATION

Air Displacement Pipets				US List Price
Model #	Description			Price
T705	5 ml Sampling System with stand			\$91.00
T710	10 ml Sampling System with stand			\$91.00
T705B	5 ml Sampling System only			\$73.00
T710B	10 ml Sampling System only			\$73.00

  

ACCESSORIES				
Polypropylene Disposable Tips				
Model #	Size ml	Qty/Pkg		US List Price
D7055	5.0	50		\$22.80
D7051	5.0	100		\$42.60
D7105	10.0	50		\$33.20
D7101	10.0	100		\$46.70