



## FI-Stream® Storage Reservoirs



Model A56290-857 with Storage Reservoir

### Complete Your Compact Distilled Water System!

- 24 hour automatic operation
- Wall or bench mounted
- Vent filter prevents air-borne contamination in reservoir

#### Product Description

#### **Reservoirs for use with FI-Stream III 4 and 8-liter models only**

##### **A56290-857: 30 Liter Capacity**

- Allows for 24 hour automatic operation.
- Vent filter prevents airborne contamination.
- Built-in overflow tubing directs water to drain in case of overflow.
- Contains special ports to feed deionization systems.
- Contains a sight tube for easy observation of water level.
- Easy-to-clean storage tank.
- Bench or wall mounting options for flexibility.

##### **A1058: 50 Liter Capacity**

- Portable unit easily moves to other lab areas needing distilled water.
- Non-automatic operation.

#### **Reservoir for use with FI-Stream II 2S models only**

##### **A1052: 8 Liter Capacity**

- Includes wall mounting bracket for easy installation.
- Portable design is ideal for transporting distilled water to different lab sites.



Model A1058 50-Liter

Model A1052 8-Liter

### ORDERING INFORMATION

Reservoirs for use with FI-Stream III 4 and 8 liter models only							
Model #	Capacity	Overall Dimensions Inches (cm)			Operating Weight Lb. (kg)	Shipping Weight Lb. (kg)	US List Price
		W	H	D			
A56290-857	30 liters	18 (46)	29.5 (75)	14.75 (37)	130 (59)	61.5 (28)	\$1,325.00
A1058	50 liters	14 (35.6)	32 (81)	14 (35)	121 (55)	11 (4.9)	\$397.00
Reservoir for use with FI-Stream 2S models only							
A1052	8 liters	8 (20)	6.5 (17)	17 (43)	7 (3.2)	25 (11)	\$346.00