



Barnstead Distillation



Classic Still Accessories



Steam still with automatic controls on floorstand connected to storage tank.



Pretreatment option for 1 and 2 gallon electric still with automatic controls mounted on wall bracket.

Fully Automatic Controls

- Provides for unattended operation.
- Stops still when storage tank is full.
- Starts still when storage tank water level falls below preset level.
- Drains still boiler at your selected 2, 4, 8, or 16 hr. intervals to reduce scale formation.
- The Controls include: control box, storage tank monitor, Inlet piping includes: solenoid valve, pressure gauge, manual regulating valve. Boiler drain includes: piping, solenoid valve, manual valve.

Wall Bracket for Classic Stills G3250, G3251, G3252, G3256, G3275

- Available for all classic stills, 1 to 10 GPH.
- The wall brackets are made of aluminum casting (G3275 made of steel).

Pretreatment Kit - 1 and 2 Gallon Still

- Uses deionization resins and carbon to eliminate cleaning and increase purity.
- Comes complete with flow meter to maintain correct flow to still.

ORDERING INFORMATION


Fully Automatic Controls Model #	Voltage	Where used	US List Price
G2100	120V 50/60 hz	1 & 2 gph 1 phase electric stills	\$2,896.00
G2101	120V 50/60 hz	2 gph 3 phase electric stills	\$2,896.00
G2110	120V 50/60 hz	5 gph electric stills, 1 & 3 phase	\$3,267.00
G2120	120V 50/60 hz	10 gph 1 phase electric stills	\$3,710.00
G2125	120V 50/60 hz	10 gph 3 phase electric stills	\$3,710.00
G2010	120V 50/60 hz	5 gph steam stills	\$4,434.00
G2020	120V 50/60 hz	10 gph steam stills	\$4,517.00

Fully Automatic Control is connected to the 120/240V Low Water Cutoff control box. 240V Low Water Cutoff is fitted with a stepdown transformer to 120V. 120V Fully Automatic Control may be connected to any Low Water Cutoff.

Wall Brackets Model #	Description	US List Price
G3250	Wall Bracket for all 1 GPH Electric Stills	\$235.00
G3251	Wall Bracket for 2 GPH Electric Still	\$235.00
G3252	Wall Bracket for 5 GPH Electric Still	\$414.00
G3256	Wall Bracket for 5 GPH Steam Still	\$414.00
G3275	Wall Bracket for 10 GPH Electric and Steam Stills	\$570.00

Pretreatment Kit Model #	Description	US List Price
G3628	Pretreatment Kit, 1 and 2 gallon still includes one cartridge	\$620.00
D8921	Replacement Pretreatment Cartridge	\$89.00

Barnstead Distillation



Puromatic Diverter for Classic Stills G3330

- Allows you to set minimum distilled water resistivity.
- Diverts substandard distilled water to drain.
- All piping constructed of inert pure tin.
- Automatic temperature compensated purity meter.

Automatic Feedwater Shut-off Valve for Classic Stills G2910

- Must be used with float feeder and fully automatic control on 5 & 10 GPH classic stills using pretreated feedwater.
- Conservation of water and automatic operation.
- To be used with float feeder (G3721) to start or stop feedwater.

Floorstand for Classic Stills G1000, G1001

- Available for all classic stills, 1 to 10 GPH.
- The floorstands are crafted of sheet metal and are provided with adjustable feet.

Float Feeder for Classic Stills

- Conserves pretreated water.
- Float mechanism controls the flow of pretreated water to boiler.
- Purity is maintained with pure tin surfaces.
- Used only on 5 and 10 gph classic stills.
- Must be used with automatic feedwater solenoid valve and fully automated controls on 5 & 10 GPH.

ORDERING INFORMATION

Puromatic Diverter				Classic Stills Automatic Shut-off Valve	
Model #	Description			Model #	Description
			US List Price		
G3330	Puromatic Diverter for Classic Stills		\$3,346.00	G2910	Automatic Feed water shut-off valve for Classic Still
					US List Price
					\$625.00
Floorstand for Classic Stills G1000, G1001					
Model #	Description	Dimensions Inches (cm) (W x D x H)	Shipping Weight Lbs. (kg)	US List Price	
G1000	Floorstand for 5 and 10 GPH Electric and Steam Stills	19.75 x 9.5 x 35 (50 x 24 x 89)	40 (18.1)	\$491.00	
G1001	Floorstand for 1 and 2 GPH Electric Stills	10.38 x 6 x 36 (26.4 x 15.2 x 91.4)	25 (11.3)	\$460.00	
Float Feeder for Classic Stills					
Model #	Description			US List Price	
G3721	Float Feeder for Classic Stills			\$1,259.00	
G2910	Automatic Feedwater Solenoid Valve			\$625.00	