



## Classic Stills Accessories



Tin-Coated Cylindrical Storage Tanks



Ventgard Air Filter

### Store Your Water and Keep It Pure!

- Store up to 200 gallons of distilled water
- Maintain purity with air filter and water seal and UV lamp

- Connections are provided for: Water inlet, distribution outlet, tank drain (not on 10 and 25 gallon)
- Connections are also provided for optional accessories: Pump, low water cut-off, pump by-pass return, still level control, ultraviolet lamp, Ventgard and water seal. Connections not used are provided with removable plugs.

#### Rectangular Tanks

- Includes ultraviolet lamp assembly, Ventgard and water seal, faucet and sight glass.
- Space-saving design.
- Wall mounted with two convenient sizes.

#### Ventgard Air Filter and Water Seal

- Ventgard Air Filter protects the stored distilled water against airborne particulates, organics and CO<sub>2</sub>.
- Water Seal works as a one-way valve to allow distilled water into storage tank. It prevents air from entering the storage tanks via the condenser atmospheric vent.

### Product Description

- Tin-lined tanks are the ideal complement to the Classic Still.
- Range in storage capacity from 10-200 gallons.
- Crafted from copper and hand-wiped with pure tin to protect the purity of the distilled water.

### Cylindrical Storage Tanks

- All tanks include a removable cover, water level sight glass, and a tin coated draw off faucet with serrated hose nipple.

### ORDERING INFORMATION

Tin-Coated Storage Tanks										
Model #	Description	Size (Gallon)	W	Dimensions Inches (cm)			Distribution Outlet NPT	Drain NPT	Shipping Weight Lb. (kg)	US List Price
				H	D					
B3027	Rectangular	25	37 (93.9)	21 (53.3)	31 (78.7)	1/2	—	150 (68)	\$6,474.00	
B3028	Rectangular	50	37 (93.9)	26 (66.0)	36 (91.4)	3/4	—	200 (90.7)	\$8,612.00	
B3043	Cylindrical	10	10.25 (26)	29 (74)	17 (43)	1/2	—	25 (11.3)	\$2,587.00	
B3045	Cylindrical	25	18 (45.7)	29 (73.7)	25 (63.5)	1/2	—	45 (20.4)	\$2,957.00	
B3046	Cylindrical	50	21.7 (55.3)	36 (91.4)	28 (71.1)	1	1/2	195 (88.5)	\$3,709.00	
B3047	Cylindrical	100	28.5 (72.4)	43.1 (109.5)	35 (88.9)	1	1/2	420 (190.5)	\$5,445.00	
B3049	Cylindrical	200	37.9 (96.3)	50.25 (127.6)	45 (114.3)	1 1/2	1	400 (181.4)	\$8,573.00	
Ventgard Air Filter										
Model #	Description									US List Price
H3111	Ventgard and 1 to 10 gph water seal									\$438.00
H3120	Ventgard air filter									\$151.00
H3130	Water seal 1 to 10 gph									\$346.00
H3131	Water seal 15 to 30 gph									\$287.00
25001	Replacement filter element									\$52.00