



Cabinetized Stills

Your Complete Distilled Water Station!

- Completely cabinetized distillation system
- Choice of two sizes: 2 or 5 GPH
- Fully automatic operation
- Automatically drains the boiling chamber every time the unit is off and every 4 hours of operation, helping to keep your still clean.

Model 210 Still

- Safety feature: interrupts power to heating elements if feed water is off.
- Optional bench stand for bench mounting of still.
- 2 GPH still and 10 gallon storage tank for the production of distilled water.

Model 525 Still

- Purity meter alerts you of distilled water purity in storage tank.
- Q-Baffle® allows for pyrogen removal.
- Optional recirculation pump and base mounting stand available.
- Storage tank is equipped with a UV lamp for bacterial control.
- 5 GPH still and 25 gallon storage tank for production of distilled water.



Model 525 Cabinetized Still



Model 210 Cabinetized Still (cover off)

Product Description

- Completely cabinetized: Includes a still and an appropriately matched tin-coated storage tank in a neat, compact and clean, unit.
- Space saving design takes up less space than conventional units, while allowing for point of use installation.
- Fully Automatic: Hands-free operation; turns still off when the tank is full, restarts still when tank empties.

PRODUCT SPECIFICATIONS

Model 525 Still						
Model #	Water Inlet (CW)	Drain Connector (WA)	W	Dimensions Inches (cm) H	D	Shipping Weight Lb. (kg)
A1085-B	Sweat Joint for 3/8" Copper tubing	Sweat Joint for 1 1/4" Copper tubing	36.5 (91)	48 (122)	18.75 (48)	425 (193)
A1085-C	Sweat Joint for 3/8" Copper tubing	Sweat Joint for 1 1/4" Copper tubing	36.5 (91)	48 (122)	18.75 (48)	425 (193)
A1085-D	Sweat Joint for 3/8" Copper tubing	Sweat Joint for 1 1/4" Copper tubing	36.5 (91)	48 (122)	18.75 (48)	425 (193)
Model 210 Still						
A1065-B	Sweat Joint for 1/2" Copper tubing	Sweat Joint for 1 1/4" Copper tubing	24 (60)	41 (103)	15.25 (39)	165 (75)
A1065-C	Sweat Joint for 1/2" Copper tubing	Sweat Joint for 1 1/4" Copper tubing	24 (60)	41 (103)	15.25 (39)	165 (75)
A1065-D	Sweat Joint for 1/2" Copper tubing	Sweat Joint for 1 1/4" Copper tubing	24 (60)	41 (103)	15.25 (39)	165 (75)

ORDERING INFORMATION

Model 525 Still							Model 210 Still						
Model #	Volts	Wire	Phase	Amps	Feedwater	US List Price	Model #	Volts	Wire	Phase	Amps	Feedwater	US List Price
A1085-B	240	2	1	55	Tap	\$20,431.00	A1065-B	230	2	1	26	Tap	\$12,968.00
A1085-C	208	4	3	36	Tap	\$20,431.00	A1065-C	208	4	3	17	Tap	\$12,968.00
A1085-D	240	3	3	33	Tap	\$20,431.00	A1065-D	240	3	3	15	Tap	\$12,968.00
							A1066	Optional Floor Stand				\$306.00	

ACCESSORIES

Model #	Description	US List Price
A1521	Floorstand 36"x18.75"x35" (92x48x89 cm)	\$576.00
A1522	Recirculating Pump (Includes A1521 floorstand)	\$3,834.00

Controls for all models require separate 120 volts.