



Analog Modular Block



Type DB17600 and DB28100 Dri-Baths

Standard Temperature Models!

- Choose from temperatures to 110°C, 115°C or 140°C
- Customize your bath with a removable block
- Available in single or three block models



Product Description

- Dri-baths can be customized by choosing the modular blocks with well sizes to fit specific applications. Modular blocks sold separately.
- Aluminum modular blocks provide excellent heat transfer to tube contents to enhance your testing results.
- All blocks include thermometer well for easily measuring block temperature.
- Must use Thermolyne modular block.

- Adjustable temperature range with superior temperature control:
 - Type DB16500 models: $\pm 0.1^\circ\text{C}$ to 110°C.
 - Type DB17600 models: $\pm 0.5^\circ\text{C}$ to 115°C.
 - Type DB28100 models: $\pm 0.5^\circ\text{C}$ to 140°C.
- Compact design saves valuable counter space.

Operation

- Two separate thermostats on the Type DB28100 provide low range or high range operation. The low range knob controls temperature. Settings from 30°C to 80°C. The high range knob controls temperature settings from 70°C to 140°C.

PRODUCT SPECIFICATIONS

Type	Description	Overall Dimensions Inches (cm)			Control Stability	Shipping Weight Lb. (kg)
		W	H	D		
DB17600	Single block	8.0 (20.3)	4.0 (10.1)	9.0 (22.8)	$\pm 0.5^\circ\text{C}$	3.5 (1.5)
DB16500	Triple block	14.5 (36.8)	4.0 (10.1)	9.0 (22.8)	$\pm 0.1^\circ\text{C}$	9.0 (4.1)
DB28100	Triple block	12.0 (30.5)	4.0 (10.1)	9.0 (22.8)	$\pm 0.5^\circ\text{C}$	7.0 (3.1)

ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hertz) Amps	Watts	Standard Calibration Range	US List Price
DB17615	120	0.6	75	Amb+5 - 115°C	\$258.00
DB17610	240	0.3	75	Amb+5 - 115°C	\$258.00
DB17610-26	220	0.3	75	Amb+5 - 115°C	\$258.00
DB16525 ¹	120	3.75	450	25°C - 110°C	\$804.00
DB16520	240	1.88	450	25°C - 110°C	\$804.00
DB28125	120	2.1	250	Amb+5 -140°C	\$544.00
DB28120	240	1.0	250	Amb+5 -140°C	\$544.00
DB28120-26	220	1.1	250	Amb+5 -140°C	\$544.00

¹UL Listed