



Modular Block Accessories



Model 2069, 2071
and 2072

- Excellent heat transfer
- Anodized aluminum provides greater heat transfer
- Choose from many varieties

Product Description

- Close contact of tubes to block walls provides excellent heat transfer.
- Blocks include thermometer well for monitoring block temperatures.
- Solid black anodized aluminum blocks for greater heat retention.
- Block puller included with all models.
- For use with all Lab-Line brand modular dri-baths.
- Single block size: 3.75" W x 2" H x 2.9" D (9.5 x 5.1 x 7.4 cm).



Model 2064 and 2065

ORDERING INFORMATION

Model #	Description	Opening Diameter Inches (cm)	US List Price
2058	Modular block: holds 30, 0.5 ml tubes and 19 .2 ml tubes	0.29 (.74) taper, 0.24 (.61) taper	\$119.00
2059	Modular block: holds 79, 0.2 ml tubes	0.24 (.61) taper	\$151.00
2064	Modular block: holds 1 titer plate, direct contact base, size of two blocks	-	\$260.00
2065	Modular block: holds 1 titer plate, size of two blocks	-	\$144.00
2066	Modular block: holds 12, 10 mm cuvettes	0.5 (1.27) square	\$100.00
2068	Modular block: holds 30, 8 mm microcentrifuge tubes (tapered for 0.5 ml centrifuge tubes)		\$84.00
2069	Modular block: holds 20, 10 mm microcentrifuge tubes (tapered for 1.5 ml microcentrifuge tubes)		\$84.00
2070	Modular block: holds 30, 6 mm test tubes	0.31 (.79)	\$79.00
2071	Modular block: holds 24, 10 mm test tubes	0.4 (1.02)	\$79.00
2072*	Modular block: holds 22, 12 mm and 13 mm test tubes	0.53 (1.35)	\$79.00
2073	Modular block: holds 12, 15 mm and 16 mm test tubes	0.68 (1.73)	\$79.00
2074	Modular block: holds eight, 20 mm test tubes	0.81 (2.06)	\$79.00
2075	Modular block: holds six, 25 mm test tubes	1.03 (2.62)	\$79.00
2076	Modular block: holds three, 25 mm test tubes; five, 12 or 13 mm test tubes and six, 6 mm test tubes	1.03 (2.62) 0.53 (1.35) 0.31 (.79)	\$79.00
2078*	Modular block: solid, use as a hot plate		\$79.00
2080	Dri-bath chamber: size of four blocks		\$213.00
2081	Modular block: holds 12, 17 mm and 18 mm test tubes	0.75 (19)	\$87.00
2083	Modular block: holds one titer plate (or 96, 0.2 ml tubes) and six, 0.5 ml tubes. size of 2 blocks	0.29 (7.5) taper	\$575.00
2085*	Stainless steel cover for models 2064 and 2065		\$33.00

*Used only for Lab-Line.