

Lab-Line Environmental Chambers



Educational Animal Study



Model 849-2

Animals in the Classroom!

- Offers full observation of animals
- Temperature from slightly above ambient to 37°C
- Hi-Lo limit thermostats for safety

Product Description

- Hardwood maple cabinet exterior.
- Easy-to-clean white enamel interior.
- Two sliding, low light transmission glass doors with key lock.
- Three adjustable wire shelves, two stainless steel pull-out trays and removable mesh bottom.

Application

- Designed to house small laboratory animals such as mice, rats or guinea pigs.

PRODUCT SPECIFICATIONS

Model	Capacity cu. ft. (l)	Shelf Area sq. ft. (sq. m)	Interior Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
			W	H	D	W	H	D	
849-2	31 (878)	22.2 (2.1)	52 (132)	43 (109)	24 (61)	56 (142)	74 (188)	28.5 (72)	800 (363)

ORDERING INFORMATION

Model	Volts	Electrical		US List Price
		Hz	Watts	
849-2	120	60	1400	\$12,648.00
849-3 ¹	120	60	1400	\$13,083.00
849-4 ²	120	60	1400	\$13,117.00

¹ Same as 849-2 except with charcoal and dust filters for indoor venting.

² Same as 849-3 except with casters.

Models 849-3 and 849-4 are designed to exhaust into the room; model 849-2 is designed for external venting (customer supplied.)