

Lab-Line
**Environmental
Chambers**



Educational Work Station



Model 890

Instant Laboratory!

- Temperature range from -20°C to 60°C
- Includes two incubators, freezer, refrigerator in addition to a work station
- Complete laboratory work station for the classroom

Product Description

- Complete classroom laboratory work station.
- Standard equipment includes: Two incubators with six perforated stainless steel shelves, Freezer with four fiberglass trays, Refrigerator with five plastic-coated wire shelves, Work center with two stainless steel shelves, chrome hosecock, sink with brass trap, 40 watt fluorescent light with switch, duplex electrical outlet.
- Cabinet exterior has maple hardwood frame.
- All cabinet interiors and doors are stainless steel.

Applications

- Tissue or cell culture incubation and preparation.
- Culture storage in the dormant state.

PRODUCT SPECIFICATIONS

Model 890	Chamber Volume cu. ft. (l)	Shelves and/or Trays	Shelf Area sq. ft. (sq. m)	External Dimensions Inches (cm)		
				W	H	D
Complete Biological Work Station	-	-	-	62 (157)	82 (208)	27 (69)
Incubator-Viewing Window, Electrical Outlet and Light	8 (227)	3	11.5 (1.1)	24 (61)	24 (61)	24 (61)
Incubator-Solid Door, Electrical Outlet and Light	8 (227)	3	11.5 (1.1)	24 (61)	24 (61)	24 (61)
Freezer	6 (173)	2	4 (0.4)	24 (61)	34.25 (87)	23.5 (60)
Refrigerator	6 (173)	2	6.25 (0.4)	24 (61)	34.25 (87)	23.5 (60)
Work Center	-	2	10 (0.9)	60 (152)	-	24 (61)

ORDERING INFORMATION

	Operating Temp Range °C	US List Price
Model 890	—	\$21,007.00
Includes:		
Incubator-Viewing Window, Electrical Outlet and Light	Slightly Above Ambient to 50°C	
Incubator-Solid Door, Electrical Outlet and Light	Slightly Above Ambient to 60°C	
Freezer	-21° to -9°C	
Refrigerator	-17° to 9°C	
Work Center	—	