

## Lab-Line Environmental Chambers



## Educational Plant



Model 846

### For Your Countertop!

- Ideal for growing plants in the classroom
- Automatic 24 hour programming
- Sliding glass doors for easy access and viewing

### Product Description

- Lighting consists of 4 fluorescent and 3 incandescent lamps for various lighting levels (customer supplied.)
- 2 interior electrical outlets to easily accommodate other electrical equipment.
- Main power switch, lamp switch, heater output, fan, thermostat and 24 timer/programmer switches all have on/off indicator lights.
- Front control panel can be removed for maintenance without disturbing experiments in progress.

- Model 846 series combine large viewing area, simplified temperature controls and an economical price.

### Applications

- Designed to accommodate plant, animal and culture studies in the classroom setting.

### PRODUCT SPECIFICATIONS

Model #	Chamber Volume cu. ft. (l)	Shelf Area sq. ft. (sq. m)	Chamber Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
			W	H	D	W	H	D	
846	14 (396)	2.2 (.20)	50 (127)	29 (74)	17 (43)	52 (132)	38 (97)	18.5 (47)	236 (107)
846-3	14 (396)	2.2 (.20)	50 (127)	29 (74)	17 (43)	52 (132)	38 (97)	18.5 (47)	236 (107)

### ORDERING INFORMATION

Model	Volts	Electrical Hz	Watts	US List Price
846	120	60	800	\$3,055.00
846-3	220	50	800	\$3,055.00

### ACCESSORIES

Model #	Description	US List Price
846-8	Chamber Shelf 6.5" x 49.5" (16.5 x 125.7 cm)	\$159.00
846-12	Mobile Chamber Stand 53" W x 32" H x 23.5" D (134.6 x 81.3 x 59.7 cm)	\$708.00
846-13	Planting Tray (Pkg. of 4)	\$126.00