



## Final Filters and Ultrafilters



Ultrafilter



D3750

- Absolute filtration down to 0.2 micron
- Pyrogen-free (Ultrafilter)
- Large surface area—long life

### Product Description

#### Final Filter Models D3750, D3751 and FL703X2

- Purity is maintained with the use of hydrophilic material eliminating recontamination.
- 0.2 micron pore size, bubble point tested to ensure absolute filtration.
- Provides bacterial free water.
- Irradiated to ensure sterility.
- Extremely low extractables.
- Autoclavable.
- Used as final filter for your most critical applications.
- Automatic air venting.
- Demonstrated log reduction value (LRV) of 9.38 when challenged with *P. Diminuta*.

#### Ultrafilter Models D4721, D4722, FL582x1 and FL550x1

- Provide pyrogen-free water for cell and tissue culture.
- Spiral-wound 10,000 dalton (Molecular weight) cutoff.
- Spiral-wound design allows for easy and complete sanitization and minimizes rinse-up time.
- Short rinse-up to 18 megohm-cm.
- Polypropylene housing.
- Kit includes fittings to connect to Analytical Grade NANOpure® (series 500 only).
- No pressure gauge or special fittings needed <0.001 eu/ml and 6+ log pyrogen removal.
- Easy to clean and replace.

#### Replacement Ultrafilter for NANOpure Diamond and EASYpure II Ultrafilter Models FL1192X1

- Hollow fiber material allows for fast rinse up to purity.
- Provides pyrogen-free water.
- Maximum flow is maintained with very little pressure drop.
- <0.001 eu/ml and 6+ Log reduction of pyrogens.

### ORDERING INFORMATION FOR FINAL FILTERS

Model #	Description	Max. Flow Rate LPM	Inlet Connection	Outlet Connection	Filters Per Pack	Max. Inlet Pressure	US List Price
D3750	Final Filter - E-pure™, NANOpure II, Dlamond, EASYpure® II. Remote Dispenser - Dlamond (D11981)	2	1/4" NPT	1/4 x 3/8" tapered hose barb	1	50 psig	\$75.00
FL703X2	EASYpure, Remote Dispenser (D8952) Final Filter-Infinity and Quad Style NANOpure	1.5	Luer connector	1/4" hose barb	1	50 psig	\$56.00
D3751	Remote Dispenser, Cross Flow Final Filter- Quad Style, NANOpure	1.5	N/A	N/A	2	50 psig	\$209.00

### ORDERING INFORMATION FOR ULTRAFILTERS

Model #	Description	Max. Inlet Pressure	Max. Flow Rate LPM	Inlet Connection	US List Price
D4721	Ultrafilter conversion kit for D4700 NANOpure	90 psig	1.5	1/2" NPFT	\$878.00
FL1192X1	Replacement Ultrafilter for NANOpure Dlamond and EASYpure II				\$340.80
FL582X1	Ultrafilter conversion kit for E-pure and NANOpure II	90 psig	1.5	1/2" NPFT	\$1,266.00
D4722	Replacement Ultrafilter Element for D4721, FL582X1 and NANOpure Bioresearch Grade	90 psig	1.5	1/2" NPFT	\$715.00
FL550X1	Replacement Ultrafilter for NANOpure Infinity and EASYpure	90 psig	1.5	1/2" NPFT	\$1,026.00