

Electrothermal Heating Mantles



Three Flask Sizes In One Mantle (CMUT)



One Mantle—Three Sizes!

- Aluminum case is chemically resistant and easy to clean
- External control allows for fume hood use
- Three-in-one design accommodates 250, 500 or 1000 ml round bottom flasks



Model CMUT1000/CEX1

Product Description

- Unique patented air flow keeps the heater cartridge cool to the touch.
- Ceramic fiber thermal insulation permits faster heater up while reducing energy consumption.
- Safety ground screen over heating elements protects the operator from shock hazard due to spills or flask breakage.
- One mantle can accommodate 250, 500 or 1000 ml round bottom flask sizes.
- Element temperature range of ambient to 450°C.
- Dual circuit switch control the lower circuit only or the upper and lower circuit together.
- Solid state external power control provides precise manual control of the output range from 3 to 100%.
- Control panel (on the controller) contains two indicator lights: "Power On" and "Heater On."
- An integral ring-stand holder on all models accommodates support stand and clamp.

ORDERING INFORMATION

Model #	Description	Voltage	US List Price
CMUT1000/CEX1 ¹	Mantle with Control	115	\$279.00
CMUT1000/CE ¹	Mantle with Control	230	\$279.00
RECMUTX1	Mantle only	115	\$112.30
RECMUT	Mantle only	230	\$119.10
MC239X1 ¹	Control only	115	\$134.00
MC239 ¹	Control only	230	\$134.00

¹ CE