

## Lab-Line Heating Mantles



# Multi-Unit Extraction Heater



## Multi-Position Square Top Extraction Heater!

- Temperature range from 40°C to 300°C
- Six individual heat controls with pilot lights
- Built-in temperature regulator
- ☆ 24 Month Warranty

### Product Description

- Six individual solid state heat controllers with pilot lights.
- Six 4.5" (11 cm) square cast aluminum hot plates, each with 300 watt heater.
- Inclined control panel mounted on stainless steel housing for ease of viewing.
- Four frame clamps, two 30" (76 cm) frame rods, two 36" (91 cm) frame rods, and twelve extraction clamps with clamp holders are included.
- A variety of optional clamps for attaching glassware to the support frame are available.
- Aluminum hot plates for rapid rise times.
- Separate indicator lights for power and heat.

### Operation

- Temperature range from 40°C to 300°C.
- Allows you to run different projects concurrently or to make simultaneous multiple runs of the same project.
- Ideal for extraction, distillation, and digestion.
- Can operate in temperatures from 15°C to 40°C.
- Glassware and power cord not included.

### PRODUCT SPECIFICATIONS

Model #	Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	
5000	30 (76)	9 (23)	11 (28)	40 (18)
5000-1	30 (76)	9 (23)	11 (28)	40 (18)

### ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)	Operating Temp Range °C (°F)	US List Price
		Watts		
5000	120	1800	40-300 (104-572)	\$2,796.00
5000-1	240	1800	40-300 (104-572)	\$2,796.00