


APPLICATIONS



Cytology



Microbiology

Histology

PRODUCTS



# HISTOLOGY

## SELECTION GUIDE

### SLIDE WARMERS

	Type	Operating Temp Range °C	Overall Dimensions Inches (cm)	Page
			W H* D	
<b>Our Large Size Economical Slide-Warmer!</b>	26000	Ambient to 65°C	14 (36)	124
			25 (64)	
			25 (64)	
<b>Higher Temperature Slide-Warmer With Drying Bars!</b>	MH6600	Ambient to 100°C	15.6 (39.6) 7 (17.8) 4 (10.2)	125
<b>Our Compact Economical Slide-Warmer!</b>	LL16100	45° to 50°C	13 (33) 4 (10.1) 4 (10.1)	126

\*Height with lid open

### FLOAT BATH

	Type	Operating Temp Range °C	Chamber Capacity Liters	Page
<b>Low-Profile Economical Float Bath!</b>	26100	Ambient to 70°C	2.3	127
<b>Top of the Line Float Bath!</b>	MH8500	Ambient to 75°C 40° to 80°C	2.5	128

### WAX DISPENSER

	Type	Operating Temp Range °C	Chamber Capacity Liters	Page
<b>The Hottest Wax Dispenser in the Market!</b>	MH8523	Ambient to 75°C	45	129

CELEBRATING  
**125** YEARS  
of INNOVATION, QUALITY and  
SERVICE SINCE 1878

800-553-0039 • 563-556-2241 • www.barnstead.com