



## Paraffin Section Mounting Bath



Model MH8515

### Top of the Line Float Bath!

- Comfortable and easy manipulation of sections
- Constant clear visibility
- Coated aluminum housing



### Product Description

- Easy manipulation of the section avoids problems with wrinkles, folds and distortions.
- 2.5 liter one piece black interior is easy to clean and provides constant clear visibility of the sample.
- Aluminum exterior with polypropylene lid is durable and easy to clean.
- Easily replaceable heating element.
- Built-in controller regulates the temperature range of 40°C to 75°C.

### Operation

- Front mounted on/off switch and temperature control knob allows for complete control of the unit.
- Clear light indicates "Power On." Red light warns of over-temperature. Amber light indicates "Heater On."
- Excess or damaged surfaces can be easily removed by skimming a filter paper over the water surface.

### Safety

- Built-in safety device shuts off the energy supply to the heater in the event of an overheat condition such as boiling dry.

### PRODUCT SPECIFICATIONS

Model #	Chamber Capacity (l)	Chamber Dimensions Inches (cm)		Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)
		H	D	H	D	
MH8515X1	2.25	3.0 (7)	8.7 (22)	5 (13)	13 (33)	6 (2.7)
MH8515	2.25	3.0 (7)	8.7 (22)	5 (13)	13 (33)	6 (2.7)

### ORDERING INFORMATION

Model #	Volts	Electrical	Watts	US List Price
		(50/60 Hertz)		
MH8515X1 <sup>1</sup>	115		175	\$1,037.00
MH8515 <sup>1</sup>	230		240	\$1,037.00

<sup>1</sup> CE