



## Paraffin Wax Dispenser



### The Hottest Wax Dispenser In the Market!

- Provides immediate "Wax on Tap"
- For use in Pathology laboratories, etc.
- Thermally insulated tap heater with independent temperature control eliminates wax solidification



Model MH8523BX1

#### Product Description

- Anodized black interior, stove enameled black and white aluminum exterior and a black metal lid are durable and easy to clean.
- Capacity of 1 gallon (4.5 liters).
- Built-in controller provides temperature control from ambient to 75°C.
- A non-drip swivel lever with replaceable filter makes the dispensing of wax easy.

#### Operation

- The clear light indicates "Power On." The red light warns of over-temperature. Amber neon lights indicate "Bath Heater On" and "Tap Heater On."
- Tap heater maintains temperature within the tap to keep the wax molten.

#### Safety

- Safety over-temperature at 105°C.
- Non-resettable thermal fuse.

#### PRODUCT SPECIFICATIONS

Model #	Chamber Capacity (l)	Chamber Dimensions Inches (cm)		Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)
		H	Diameter	H	Diameter	
MH8523BX1	(4.5)	8.25 (15)	8.25 (15)	12 (30)	11 (28)	10 (4.5)
MH8523B	(4.5)	8.25 (15)	8.25 (15)	12 (30)	11 (28)	10 (4.5)

#### ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hertz)	US List Price
		Watts	
MH8523BX1 <sup>1</sup>	115	175	\$2,318.00
MH8523B <sup>1</sup>	230	175	\$2,318.00

<sup>1</sup> CE