

Lab-Line
Hot Plate and
Hot Plate Stirrers



Pyroplate[®] Aluminum Top Hot Plates



Model 1180



Product Description

- Compact design conserves valuable bench top space.
- Aluminum top provides rapid heat transfer and excellent uniformity.
- Durable cold-rolled steel body construction.

Safe and Economical!

- Adjustable temperatures up to 382°C
- Reflective aluminum top
- ☆ 24 Month Warranty

Operation

- Heats beakers and flasks up to 1 l capacity.
- Pilot light alerts user when unit is heating.

Applications

- Ideal for general lab heating
- Sample drying
- Heating TLC plates
- Liquid evaporation

PRODUCT SPECIFICATIONS

Type	Heating Surface Inches (cm)		Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	W	H	D	
1180	8 (20)	7 (18)	8.25(21)	4.625 (12)	8.625 (22)	8 (4)

ORDERING INFORMATION

Model	Electrical (50/60 Hertz)		Maximum Operating Temp		US List Price
	Volts	Watts	°C	°F	
1180	120	500	382	720	\$372.00
1180-1CE ¹	240	500	382	720	\$372.00