

Lab-Line  
Hot Plate and  
Hot Plate Stirrers



## Aluminum Top Pyro-Magnestir® Hot Plate Stirrers



Model 1179

### Versatile Aluminum Tops!

- \* Designed for continuous use, to handle fluids of varying viscosities in up to 1 l capacity
- \* Multiple stirring option
- ☆ 24 month warranty



### Product Description

- Available in 5 models to accommodate most laboratory applications requiring mixing and temperature control.
- Aluminum tops provide uniform heat transfer to samples.
- Models 1179 and 1266 heat and stir flasks up to 1 l.
- Model 1263 simultaneously stirs and heats up to four 1 l flasks.
- Model 1268 simultaneously stirs and heats up to six 1 l flasks.
- Model 1287 provides independent stirring control for up to six 1 l flasks.

### Application

- Gentle heating

### PRODUCT SPECIFICATIONS

Model #	Platform Dimensions Inches (cm)			Overall Dimensions Inches (cm)		# of Stirring Positions	Shipping Weight Lb. (kg)
	W	D	W	H	D		
1179	7 (18)	8 (20)	8.25 (21)	5 (13)	8.625 (22)	1	11 (5)
1179-1CE	7 (18)	8 (20)	8.25 (21)	5 (13)	8.625 (22)	1	11 (5)
1263	10 (25)	10 (25)	11 (28)	7.125 (18)	10 (25)	4	23 (10)
1263-1CE	10 (25)	10 (25)	11 (28)	7.125 (18)	10 (25)	4	23 (10)
1266	6 (15)	6 (15)	7.5 (19)	4.75 (12)	6 (15)	1	7 (3)
1266-1CE	6 (15)	6 (15)	7.5 (19)	4.75 (12)	6 (15)	1	7 (3)
1268	12 (30)	18 (46)	18 (46)	7.5 (19)	14 (36)	6	42 (19)
1268-1CE	12 (30)	18 (46)	18 (46)	7.5 (19)	14 (36)	6	42 (19)
1287	12 (30)	18 (46)	18 (46)	7.5 (19)	14 (36)	6	40 (18)
1287-1CE	12 (30)	18 (46)	18 (46)	7.5 (19)	14 (36)	6	40 (18)

### ORDERING INFORMATION

Model #	Electrical (50/60 Hertz)		Operating Temp Range		Stirring Speed Range RPM	US List Price
	Volts	Watts	°C	°F		
1179 <sup>1,2</sup>	120	575	80-371	176-700	60-1650	\$500.00
1179-1CE <sup>3</sup>	240	575	80-371	176-700	60-1650	\$500.00
1263 <sup>2</sup>	120	700	80-320	176-608	100-1650	\$1,419.00
1263-1CE <sup>2,3</sup>	240	700	80-320	176-608	100-1650	\$1,419.00
1266 <sup>3</sup>	120	500	80-382	176-720	100-1650	\$409.00
1266-1CE <sup>3</sup>	240	500	80-382	176-720	100-1650	\$409.00
1268 <sup>*2</sup>	120	2200	80-320	176-608	100-1200	\$2,290.00
1268-1CE <sup>3</sup>	240	2200	80-320	176-608	100-1200	\$2,290.00
1287 <sup>2</sup>	120	2200	80-371	176-700	60-1200	\$2,400.00
1287-1CE <sup>3</sup>	240	2200	80-371	176-700	60-1200	\$2,400.00

<sup>1</sup>UL, <sup>2</sup>ETL, <sup>3</sup>CE

\* Has a 20-amp cord set.