

## Lab-Line Incubators



## Compact Low Cost



Model 150



Model 100



## Compact Low-Cost for the Clinic or Classroom!

- Aluminum chamber provides excellent temperature uniformity
- Choice of three chamber sizes
- Choose from metal or clear Plexiglas® door
- Small footprint saves valuable bench space
- ☆ 24 Month Warranty

•Plexiglas is a registered trademark of Rohm and Haas Corporation

**Includes FREE Thermometer**

### Product Description

- Corrosion-resistant aluminum chambers with 1" (2.5 cm) fiberglass insulation provide excellent temperature uniformity. (Model 150 does not include 1" insulation.)
- Plexiglas door option permits unobstructive viewing of samples and eliminates temperature changes from opening and closing the door.
- Rugged metal door with positive latch provides an excellent seal ensuring temperature stability and uniformity.

### Operation

- Bimetallic thermostat controls heat uniformity in chamber, allowing adjustment of temperature of the unit to match thermometer setting.

- Red, spirit-filled thermometer contains no mercury and is accurate to NIST tolerances.
- Pilot light indicates when heaters are energized, allowing you to be confident that the unit is operating.
- Specially designed shielded heating elements eliminate safety hazards caused by spills.

### Applications

- Bacteriological cultures
- Hematology studies
- Chick embryo studies

### PRODUCT SPECIFICATIONS

Model #	Operating Temp Range °C	Temperature Control °C	Temperature Uniformity °C	Chamber Volume cu. ft. (l)	Chamber Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
					W	H	D	W	H	D	
100, 101, 105, 106	Ambient +5° to 40°C	±1°C at 37°C	±1°C at 37°C	0.67 (19)	12 (30)	8 (20)	12 (30)	13.5 (34)	11.8 (30)	14.5 (37)	15 (7)
120, 121, 125, 126	Ambient +5° to 60°C	±1°C at 37°C	±1°C at 37°C	1.0 (28)	12 (30)	12 (30)	12 (30)	13.5 (34)	15.7 (40)	13.9 (35)	27 (12)
150, 152	Ambient +5° to 45°C	±3°C at 37°C	±3°C at 37°C	1.32 (37)	13.44 (35)	12.68 (32)	12.65 (35)	13.5 (37)	15.7 (40)	13.9 (35)	33 (15)

### CHOOSE METAL OR PLEXIGLAS DOORS

Model #	Door	Shelves (Perforated Aluminum)	Model #	Door	Shelves
100, 101	Metal	One fixed shelf	105, 106	Plexiglas	One fixed shelf
120, 121	Metal	One removable shelf	125, 126	Plexiglas	One removable shelf
			150, 152	Plexiglas	Two fixed shelves

### ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)		US List Price	Model #	Volts	Electrical (50/60 Hz)		US List Price
		Amps	Watts				Amps	Watts	
100 <sup>3,4</sup>	120	0.8	100	\$343.00	121	240	0.4	100	\$343.00
101	240	0.4	100	\$388.00	125 <sup>3,4</sup>	120	0.8	100	\$388.00
105 <sup>3,4</sup>	120	0.8	100	\$347.00	126	240	0.4	100	\$347.00
106	240	0.4	100	\$388.00	150 <sup>*3,4</sup>	120	1.3	150	\$388.00
<sup>3</sup> UL, <sup>4</sup> CUL					152	240	0.6	150	\$365.00

\*Model 150 Educational incubator for egg incubation. Unit not intended for clinical use, microbiological or biochemical studies.