




Set-Up Frame Sets



Set-Up Frame Sets

Strong and Versatile Support Systems!

- Easy construction in almost any rectangular shape and size
- Two sizes of frame set-ups
- Eight sizes of frame rods
- ☆ 24 Month Warranty

Product Description

- Frame sets come with all components necessary for a complete frame assembly.
- Extruded solid aluminum .5" (1.3 cm) diameter rods.
- Rigid frame supports allow the addition or removal of supports and apparatus without needing to dismantle the entire set-up.

Applications

- Wide assortment of rod lengths available, suitable for many laboratory applications.
- Rigid frame supports are especially useful for glassware and apparatus assemblies in organic chemistry laboratories.

ORDERING INFORMATION

Model #	Description	Dimensions, Inches (cm)		US List Price
		H	W	
7090	Includes 4 ea. 7100-1, 2 ea. 7100-3, 8 ea. 7100-4, and 22 ea. 7050-1	24 (61)	24 (61)	\$338.00
7096	Includes 6 ea. 7100-1, 6 ea. 7100-2, 2 ea. 7100-4, 7 ea. 7100-6, 5 ea. 7100-7, 49 ea. 7050-1 and 8 ea. 7051	48 (122)	72 (183)	\$768.00

ACCESSORIES

Model #	Description	Length Inches (cm)	US List Price
7100-1	Alumaloy Rods	1.64 (4.1)	\$8.00
7100-2	Alumaloy Rods	12 (30.5)	\$10.00
7100-3	Alumaloy Rods	18 (45.7)	\$11.00
7100-4	Alumaloy Rods	24 (61)	\$15.00
7100-5	Alumaloy Rods	36 (91.4)	\$17.00
7100-6	Alumaloy Rods	48 (122)	\$25.00
7100-7	Alumaloy Rods	72 (182.9)	\$31.00
7100-8	Alumaloy Rods	96 (243.8)	\$36.00
7050-1	Frame Clamp		\$12.00
7051	Frame Foot		\$17.00