

Alumaloy® Clamps

☆ 24 Month Warranty

Model 7049—Chain Clamp



Chain Clamp

- 13" chain locks over a hook on the adjustable portion of the collar.
- Accommodates beakers with outside diameter of 2.5" to 6.5".
- Stainless steel shaft is 7.25" L x 0.43" dia. (18.4 x 1.6 cm).
- Total length 11.5" (29.2 cm).

Model 7051—Frame Foot



Frame Foot

- Designed to support frame set-ups.
- Accepts .5" (1.3 cm) diameter rods.
- Three countersunk screw holes in foot.
- Allen set screws included.
- Overall diameter 2" (5 cm).

Model 7057—Double Support Base & Rod



Double Support Base and Rod

- White glazed porcelain base measures 13" L x 7" W (33 x 17.8 cm).
- Base features a drip trough around the rim to prevent lab bench spills.
- Solid aluminum rod is threaded at one end, measures 23" H x .5" dia. (58.4 x 1.3 cm).
- Model 7058 includes the Double Burette Holder (Model 7086) with the Double Support Base and Rod (Model 7057).

Model 7086—EMR Double Burette Holder



"ERM" Double Burette Holder

- Quickly and securely supports two burettes—from micro size to 100 ml.
- Built-in meniscus reader.
- Holder will not obscure graduation lines and numbers.
- Sturdy, corrosion-resistant Lexglas® construction.
- Overall dimensions 9.5" W x 5" H (24.1 x 12.7), weighs only 6 oz.

Model 7078—Water Bath Clamp



Water Bath Clamps

- Support thermometers and other similar apparatus over water baths.
- Alumaloy construction with three-prong, vinyl-covered jaws.
- Model 7078 has a large clamp with rotating jaws, 1" (2.54) jaw capacity.
- Model 7079 has a small clamp with .5" (1.3 cm) jaw capacity, jaws do not rotate.
- 2.5" (6.35 cm) from apparatus center to center of support.

ORDERING INFORMATION


Model #	Description	Jaw Capacity Inches (cm)	Rod Capacity Inches (cm)	Total Length Inches (cm)	US List Price
7049	Chain clamp	2.5 - 6.5	up to 0.5 (1.3)	11.5 (29.2)	\$34.00
7051	Frame foot	2.5 - 6.5	up to 0.5 (1.3)	11.5 (29.2)	\$17.00
7057	Double support base and rod	2.5 - 6.5	up to 0.5 (1.3)	11.5 (29.2)	\$71.00
7058	Double support base and rod with burette holder	2.5 - 6.5	up to 0.5 (1.3)	11.5 (29.2)	\$94.00
Model #	Description	Jaw Capacity Inches (cm)	Rod Capacity Inches (cm)	Total Length Inches (cm)	US List Price
7086	Double burette holder	0.5 (1.3)	up to 0.5 (1.3)	3.5 (8.9)	\$32.00
7079	Small water bath clamp	0.5 (1.3)	up to 0.5 (1.3)	3.5 (8.9)	\$19.00
7078	Large water bath clamp	1 (2.54)	up to 0.5 (1.3)	3.5 (8.9)	\$26.00




Alumaloy® Clamps



Model 7054—Duo-Clamp



Model 7064—All Angle Swivel Clamp



Model 7077—Thermometer Clamp

☆ 24 Month Warranty

Duo-Clamp

- Add attachments to support framework without disturbing previous set-up.
- Joins two rods in a "T" connection at any angle within a 360° arc.
- Alumaloy construction.
- Allen set screws and bead wrench included.
- 1.1" (2.8 cm) overall diameter.

All-Angle Swivel Clamp

- Connect rods, extension clamps, frames or other apparatus up to .5" diameter at any angle.
- Two 45° angle ribbed thumbscrews rotate on a connecting rod.
- Total length 3" (7.62 cm).

Thermometer Clamp

- Support thermometers, glass tubes or similar apparatus up to .5" diameter.
- Stainless steel clamping jaws have .5" (1.3 cm) adjustable extension and rotate 360°.
- 3.75" (9.5 cm) from thermometer to center of support; total length 5.5".

ORDERING INFORMATION

Model #	Description	Rod Capacity Inches (cm)	US List Price
7054	Duo-clamp	up to 0.5 (1.3)	\$82.00
7064	All-angle swivel clamp	up to 0.5 (1.3)	\$28.00
7077	Thermometer clamp	up to 0.5 (1.3)	\$21.00