



Laboratory Jacks and Saf-T-Shields



Model 25

When Lab Equipment Needs A Lift!

- Stable in all positions
- Available in three sizes
- Sturdy, corrosion-resistant anodized aluminum construction
- ☆ 24 Month Warranty

Multi-Jacks®

- Provides easy raising and lowering of many different types of instruments or apparatus, such as heaters, hot plates, water baths or glassware.
- Accurate, slippage-free vertical adjustment.
- Stable in all positions.
- Stepless screw-type control smoothly raises or lowers loading platform.
- Supports various types and sizes of glassware, apparatus and instrumentation at specific heights as required to accommodate experimental set-ups.

ORDERING INFORMATION

Model #	Size of Plate Inches (cm)		Adjustable From-To Inches (cm)	Max. Load At Middle Lb. (kg)	Shipping Weight Lb. (kg)	US List Price
	W	L				
24	4.7 (12)	5.5 (14)	2.4 (6) to 9.6 (24.3)	110 (50)	1.86 (0.85)	\$186.00
25	4.8 (12.2)	5.9 (15)	3.1 (8) to 10.8 (27.8)	176 (80)	2.60 (1.18)	\$222.00
26*	7.9 (20.0)	7.9 (20)	3.1 (8) to 10.9 (27.8)	176 (80)	4 (1.82)	\$304.00

*Includes two (7.9" x 7.9") platforms



Model 80

Accident Prevention!

- Impact resistant shield protects lab workers
- Weighted base prevents tip over
- For potentially dangerous chemical reactions

Saf-T-Shield

- 0.2" (.5 cm) thick Lexan® shield provides clear visibility and protection.
- Helps protect lab workers or students from the initial force of an explosion, sheet of flame, flying glass or chemical spill.
- Shield is 36" (91.4 cm) high, adjusts from 14" to 32" (35.6 to 81.2 cm) wide.
- Can fully enclose most lab set-ups such as boiling flasks, distillation or reflux apparatus—use anywhere there's a risk of a reaction getting out of hand.

Lexan® is a registered trademark of GE Plastics.

ORDERING INFORMATION

Model #	Dimensions Fully Extended Inches (cm)		Shipping Weight Lb. (kg)	US List Price
	W	H		
80	32(81.3)	36(91.4)	37(16.8)	\$680.00