



Kwik-Set® Lab-Chron® Timers



Model 1405

Easy-To-Use Timer!

- Push-button reset for ease of use
- Angled display for easy viewing
- Precision motor provides accurate timing
- ☆ 24 Month Warranty

Kwik-Set Lab-Chron Timers

- Timing is accurate—independent of voltage fluctuations or electrical line load.
- Large illuminated start-stop switch facilitates sequential timing.
- Five-digit sealed counter with large, easy-to-read numbers.
- High torque synchronous motor for instant starts.
- Electric brake for instantaneous stops without overshoot.
- Compact, lightweight and portable.
- Non-skid rubber feet.

Applications

- Ideal for timing growth processes, rates of chemical reactions, flow of liquids and gases and histological preparations.
- Can be used in viscosity testing and a wide range of other industrial applications.

PRODUCT SPECIFICATIONS

Type	Time Range (Increments)	Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
		W	H	D	
1405	0.1 to 9999.9 sec. (0.1 sec.)	7.37 (18.7)	4.37 (11.1)	4 (10.2)	5 (2.3)

ORDERING INFORMATION

Model #	Electrical		US List Price
	Volts	Hz	
1405	120	60	\$306.00