




Imperial V Controllers



Single Set Point Controller

Single Set Point Controller

- Simultaneous display of set temperature and actual measured temperature.
- Adjustable high limit safety thermostat.
- Standard controller for the Imperial V digital, mechanical ovens.



Single Set Point Controller With Optional Timer

Single Set Point Controller with Optional Timer

- Same features as single set point controller with additional 24 hour timer.
- Timer can be set for 24 hours in 15 minute increments.



Four Program Controller

Four Program Controller

- Programmable—four programs, each with eight ramps/dwells.
- Option “C”—ramp from as low as 1°C per minute, dwell from 0 to 999.9 hours.
- Option “D”—ramp from as low as 1°C per minute, dwell from 0 to 999.9 minutes
- Adjustable high limit safety thermostat.



Stacking rack

Multi-Program Controller

(Not pictured)

- Programmable—four programs with six steps each.
- Ramp adjustable from as low as 0.01°C per minute.
- Adjustable high limit safety thermostat.

Stacking Racks

- Allows two ovens the same size to be stacked and operated safely. Fan removes heat generated by bottom unit, preventing damage to top oven.

ACCESSORIES

Controller Model*	Description	US List Price
B	Multi-program controller with 4 programs with six steps.	\$2,206.26
C	Four program controller with 4 programs; each program has 8 ramps and 8 dwells. Timer allows you to ramp in minutes and dwell in hours.	\$2,671.31
D	Four program controller with 4 programs; each program has 8 ramps and 8 dwells. Timer allows you to ramp in minutes and dwell in minutes.	\$2,671.31
T	24 hour timer (not for use with B, C or D options) for use with single set point controller (120V/60Hz only).	\$129.78
H	24 hour timer (not for use with B, C or D options) for use with single set point controller (240V/50Hz only).	\$129.78
I	24 hour timer (not for use with B, C or D options) for use with single set point controller (240V/60Hz only).	\$129.78

*A model 3488M with a single program controller which ramps in minutes and dwells in minutes would be ordered as 3488MB. Controller price must be added to the oven price.

The US List price for a 3488MB : \$3,620.00 + \$2,206.26 = \$5,826.26