

Lab-Line
Refrigerators/Freezers



Flammable Material Storage



Model 3552-10



Model 3559-10



Model 3556-4

Safe Cold Storage!

- Designed for safe cold storage of volatile materials
- Engineered with no internal electrical components which could trigger an explosion
- Meets OSHA and National Fire Protection Association safety requirements



Product Description

- Provide a wide variety of units to keep your lab safe from potentially destructive explosions.
- Available in many sizes and styles, including uprights, combination refrigerator/freezers, and undercounter models.
- Thermostat located on exterior of cabinet to eliminate interaction with stored materials.
- Magnetic vinyl gaskets for positive door seal.
- CFC-free insulation and coolant.

Model 3551-10

- 10.1 cubic ft. refrigerator/freezer with separate doors.
- Three adjustable shelves and four door shelves in refrigerator section.
- Hydraulic thermostat; adjustable natural air flow vent.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3552-10

- 20 cubic ft. upright freezer.
- Built-in evaporators in shelves for even cooling.
- Four permanently positioned wire shelves.
- Enameled steel exterior with ABS plastic interior.
- Door equipped with key lock.
- Manual defrost.
- UL listed.

Model 3556-4

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed wire shelves, three door shelves.
- Aluminum evaporator coil.
- Manual defrost.
- UL listed.

Lab-Line Refrigerators/Freezers



Model 3556-6

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Hydraulic thermostat.
- Two adjustable wire shelves, three door shelves.
- Manual defrost.
- UL listed.

Model 3556

- 5.6 cubic ft. undercounter refrigerator.
- Hydraulic thermostat.
- Three adjustable wire shelves plus three door shelves.
- Aluminum evaporator coil.
- Manual defrost.
- Key locking door.
- UL listed.

Model 3559-10

- 10 cubic ft. single door refrigerator with interior freezer.
- Three adjustable shelves plus four door shelves.
- Hydraulic thermostat.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3566-10

- 20 cubic ft. one-door refrigerator.
- Five adjustable shelves.
- Enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3560-10

- 20 cubic ft. corrosion-resistant refrigerator designed to store caustic materials and strong acids.
- Stainless steel evaporators.
- Five adjustable stainless steel shelves.
- Enameled steel interior and exterior with stainless steel door liner for greater corrosion resistance.
- Door comes with built-in key lock.
- Manual defrost.
- UL listed.

PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3551-10, 3551-11	19.5 (50)	39 (99)	18 (46)	18 (46)	12 (32)	16.25 (41)	10.1	23.5 (60)	60 (152)	29.25 (74)	150 (68)
3552-10, 3552-11		N/A		27 (69)	57 (145)	18 (46)	20.6	32 (81)	71.5 (182)	37 (94)	240 (109)
3556	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23.5 (59)	33.5 (85)	28 (71)	120 (55)
3556-4		N/A		18 (46)	27.5 (70)	15 (38)	5.6	21.25 (54)	33.5 (85)	29.25 (74)	120 (55)
3556-6	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	29.25 (74)	120 (55)
3559-10, 3559-11	19.75 (50)	32 (81)	17.5 (45)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	27.5 (70)	150 (68)
3560-10, 3560-11	27 (69)	55 (140)	20 (51)		N/A		20.6	32 (81)	71.5 (182)	32.5 (83)	280 (127)
3566-10, 3566-11	27 (69)	55 (140)	20 (51)		N/A		20.6	32 (81)	71.5 (182)	32.5 (83)	280 (127)

ORDERING INFORMATION

Model #	Electrical			Temperature Range °C (°F)		US List Price
	Volts	Hz	Watts	Refrigerator Section	Freezer Section	
3551-10 ¹	120	60	230	-1° to 8° (30° to 46°)	-20° to -12° (-4° to -10°)	\$3,106.00
3551-11	240	50	230	-1° to 8° (30° to 46°)	-20° to -12° (-4° to -10°)	\$3,106.00
3552-10 ¹	120	60	516	N/A	-23° to -5° (-9° to 23°)	\$3,534.00
3552-11	240	50	516	N/A	-23° to -5° (-9° to 23°)	\$3,534.00
3556 ¹	120	60	160	2° to 7° (36° to 45°)	N/A	\$2,646.00
3556-4 ¹	120	60	200	N/A	-26° to -15° (-15° to 5°)	\$2,646.00
3556-6 ¹	120	60	160	-9° to 7° (15° to 44°)	-15° to 1° (5° to 33°)	\$2,646.00
3559-10 ¹	120	60	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)	\$2,893.00
3559-11	240	50	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)	\$2,893.00
3560-10 ¹	120	60	480	-5° to 15° (23° to 59°)	N/A	\$5,583.00
3560-11	240	50	480	-5° to 15° (23° to 59°)	N/A	\$5,583.00
3566-10 ¹	120	60	480	-5° to 15° (23° to 59°)	N/A	\$3,164.00
3566-11	240	50	480	-5° to 15° (23° to 59°)	N/A	\$3,164.00

¹ UL listed