

Lab-Line Refrigerators/Freezers



General Purpose



Model 3767



Model 3772



Model 3751

Sizes for Your Every Need!

- Available as all-refrigerator, combo refrigerator/freezer, upright freezer, and chest freezer models in many popular sizes
- Adjustable shelves or baskets included
- CFC-free insulation and coolant



Product Description

Model 3750

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Adjustable thermostatic control.
- Two adjustable wire shelves, three door shelves.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3751

- 5.6 cubic ft. undercounter refrigerator.
- Manual defrost.
- Adjustable hydraulic thermostat control.
- Three adjustable wire shelves, three door shelves.
- Key lock.
- UL listed.

Model 3752

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed shelves.
- Manual defrost.

Model 3753

- 4 cubic ft. undercounter freezer.
- Three adjustable shelves.
- External controls include an adjustable thermostat with deep freeze setting and an over temperature warning with visual alarm.
- Built-in keylock.
- Suitable for -25°C storage requirements.
- Manual defrost.
- UL listed.

Lab-Line Refrigerators/Freezers



Model 3761

- 7 cubic ft. chest-style freezer.
- Includes storage basket.
- Textured enameled steel exterior.
- Key lock in lid.
- Manual defrost.

Model 3762

- 18 cubic ft. chest-style freezer.
- Largest chest-style model.
- Storage basket included.
- Textured enameled steel exterior.
- Key lock in lid.

Model 3763

- 14 cubic ft. upright two-door combination refrigerator/freezer.
- Hydraulic thermostat control with separate “cold air flow” control for freezer.
- Two adjustable shelves in refrigerator.
- Magnetic vinyl gaskets for positive door seal.
- Automatic defrost.
- UL listed.

Model 3764

- 10 cubic ft. upright combination refrigerator/freezer.
- Single-door model with internal freezer compartment.
- Three adjustable wire shelves.
- Manual defrost.

Model 3766

- 20 cubic ft. frost-free refrigerator.
- Five permanently mounted shelves plus storage basket.
- Hydraulic thermostat.
- Manual defrost.
- Locking door comes standard.

Model 3767

- 20 cubic ft. upright freezer.
- Four fixed wire shelves and six door shelves.
- Door comes with built-in key lock.
- Manual defrost.
- Enameled steel interior and exterior.

Model 3770

- 1.9 cubic ft. bench top combination refrigerator/freezer.
- Single door with interior freezer compartment.
- One removable shelf included.
- Manual defrost.
- UL listed.

Model 3771

- 8.8 cubic ft. combination refrigerator/freezer with separate doors.
- Four adjustable shelves in refrigerator and 2 in freezer section.
- Automatic defrost.
- Both doors include key lock for security.
- UL listed.

Model 3772

- 11.8 cubic ft. combination refrigerator/freezer with separate doors.
- Two independent compressors allow for controlling both sections separately.
- Four adjustable shelves in refrigerator and 3 pull out shelves in freezer section.
- Refrigerator: automatic defrost Freezer: manual defrost.
- Internal controls include power supply indicator, freezer thermometer, deep freeze setting for freezer and warning with visual alarm (red LED).
- UL listed.

PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions* Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3750 ^{1,2} , 3750-1	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3751 ^{1,2} , 3751-1	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23.5 (59)	33.3 (85)	24.5 (62)	120 (55)
3752 ^{1,2} , 3752-1		N/A		18 (46)	27.5 (70)	15 (38)	5.60	21.5 (55)	33.5 (85)	26 (66)	120 (55)
3753 ^{1,2}		N/A		17.5 (45)	24 (61)	16 (41)	4	24 (61)	34 (86.4)	24 (61)	100 (45.3)
3761, 3761-1		N/A		26.25 (67)	28 (71)	16.5 (42)	7	31.5 (80)	35 (89)	24.5 (62)	125 (57)
3762, 3762-1		N/A		50.25 (128)	28 (71)	22.25 (57)	19.7	61.25 (156)	35 (89)	29.5 (75)	200 (91)
3763 ^{1,2} , 3763-1	25 (64)	36.5 (93)	15.5 (39)	24 (61)	15.25 (39)	16 (41)	14.6	28 (71)	61.5 (156)	29 (74)	190 (86)
3764, 3764-1	19.75 (50)	32 (81)	17.5 (45)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	25 (64)	150 (68)
3766, 3766-1	27 (69)	57 (145)	17.75 (45)		N/A		20	32 (81)	71 (180)	29.5 (75)	280 (129)
3767, 3767-1		N/A		27 (69)	57 (145)	18 (46)	20	32 (81)	71.5 (182)	29 (74)	240 (109)
3770 ^{1,2}	16.25 (41)	15.25 (39)	11 (28)	6.25 (16)	6 (15)	8.25 (21)	1.7	19 (47)	19 (48)	19.5 (50)	52 (24)
3771 ^{1,2}	18.5 (47)	35.5 (90)	14 (36)	18.5 (47)	13 (33)	10.5 (27)	8.8	24 (61)	61 (155)	23 (59)	160 (73)
3772 ^{1,2}	19 (48)	35 (89)	16 (41)	17 (43)	28 (71)	16 (41)	11.8	23.5 (60)	78.5 (200)	25 (64)	200 (91)

Upon installation, we recommend clearance of 3" in rear, 2" on each side and 12" above unit.

¹ UL, ² CUL