

Lab-Line
Refrigerators/Freezers



General Purpose



Model 3753



Model 3763

ORDERING INFORMATION

Model #	Volts	Electrical		Temperature Range °C (°F)		US List Price
		Hz	Watts	Refrigerator Section	Freezer Section	
3750 ^{1,2}	120	60	160	-9° to 7° (15° to 44°)	-15° to 1° (5° to 33°)	\$1,334.00
3750-1	240	50	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)	\$1,334.00
3751 ^{1,2}	120	60	160	2° to 7° (36° to 45°)	N/A	\$1,334.00
3751-1	240	50	160	2° to 7° (36° to 45°)	N/A	\$1,334.00
3752 ^{1,2}	120	60	130	N/A	-26° to -15° (-15° to 5°)	\$1,334.00
3752-1	240	50	130	N/A	-26° to -15° (-15° to 5°)	\$1,334.00
3753 ^{1,2}	120	60	144	N/A	-35° to -20° (-31° to -4°)	\$1,335.00
3761	120	60	276	N/A	-21° to -4° (-5° to 24°)	\$1,267.00
3761-1	240	50	276	N/A	-21° to -4° (-5° to 24°)	\$1,267.00
3762	120	60	516	N/A	-28° to -18° (-18° to 0°)	\$2,174.00
3762-1	240	50	516	N/A	-28° to -18° (-18° to 0°)	\$2,174.00
3763 ^{1,2}	120	60	420	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)	\$1,617.00
3763-1	240	50	420	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)	\$1,617.00
3764	120	60	360	1° to 10° (34° to 194°)	-17° to -6° (-1° to 21°)	\$1,459.00
3764-1	240	50	360	1° to 10° (34° to 194°)	-17° to -6° (-1° to 21°)	\$1,459.00
3766	120	60	480	-5° to 9° (23° to 59°)	N/A	\$2,635.00
3766-1	240	50	480	-5° to 9° (23° to 59°)	N/A	\$2,635.00
3767	120	60	870	N/A	-23° to -5° (-9° to 23°)	\$1,660.00
3767-1	240	50	870	N/A	-23° to -5° (-9° to 23°)	\$1,660.00
3770 ^{1,2}	120	60	132	0° to 7° (32° to 45°)	-5° to 0° (23° to 32°)	\$435.00
3771 ^{1,2}	120	60	120	-6° to 7° (21° to 45°)	-20° to 8° (-4° to 18°)	\$1,400.00
3772 ^{1,2}	120	60	360	2° to 7° (36° to 45°)	-28° to 16° (-18° to 3°)	\$2,500.00

¹ UL, ² CUL