

Lab-Line  
Refrigerators/Freezers



## Chromatography Refrigerators



Model 3796



### Ready to Use!

- Three models available—23, 47, and 69 cu. ft.
- CFC-free insulation and coolant
- ☆ 60 Month Warranty on Compressor

### Product Description

- Temperature control of  $\pm 0.5^\circ$  at  $4^\circ\text{C}$ .
- Adjustable thermostat control from  $1^\circ$  to  $10^\circ\text{C}$ .
- Triple tempered glass doors.
- 2" foamed-in-place polyurethane insulation.
- Two 2" diameter PVC access ports for gassing applications.
- White vinyl exterior coating is durable and chemically resistant.
- Interior construction is .040 aluminum with baked-on enamel finish.
- Vapor-proof outlets.
- Manual defrost.
- Additional accessories available, contact Customer Service for ordering information.

### PRODUCT SPECIFICATIONS

Model #	Capacity Cu. Ft.	Access Ports 2" Diameter	Interior Dimensions Inches (cm)			Exterior Dimensions* Inches (cm)			Shipping Weight Lb. (kg)
			W	H	D	W	H	D	
3790, 3790-1	23	1	24 (61)	61.5 (157)	27 (69)	27.25 (70)	78.25 (199)	29.5 (75)	311 (141)
3793, 3793-1	47	2	50.75 (129)	61.25 (156)	24.5 (63)	54 (138)	78.25 (199)	29.5 (75)	495 (225)
3796, 3796-1	69	2	75 (191)	61.5 (157)	27 (69)	78 (199)	79.25 (202)	29.5 (75)	678 (308)

\*Exterior height does not include leveling feet. Add 6" for feet.

### ORDERING INFORMATION

Model #	Electrical		Duplex Outlets	Chromatography		Doors	Compressor	US List Price
	Volts	Hz		Rods	Shelves			
3790 <sup>1</sup>	120	60	1	N/A	4	1 (swinging)	1/3 hp	\$3,954.00
3790-1CE <sup>1,2</sup>	240	50	1	N/A	4	1 (swinging)	1/3 hp	\$3,954.00
3793 <sup>1</sup>	120	60	2	1	8	2 (sliding)	1/2 hp	\$6,164.00
3793-1CE <sup>1,2</sup>	240	50	2	1	8	2 (sliding)	1/2 hp	\$6,164.00
3796 <sup>1</sup>	120	60	2	2	12	3 (sliding)	1/2 hp	\$7,747.00
3796-1	240	50	2	2	12	3 (sliding)	1/2 hp	\$7,747.00

<sup>1</sup> UL, <sup>2</sup> CE

### ACCESSORIES

Model #	Description	Model #	Description	US List Price
C	Castors	3790-8	Shelf for 3790 and 3790-1CE	\$142.00
CR	Chromatography Rods	3793-8	Shelf for 3793 and 3793-1CE	\$142.00
DT	Digital Temperature Readout	3796-8	Shelf for 3796 and 3796-1CE	\$142.00
P	Access Ports			
V	Vapor-proof Duplex Outlet			

When ordering, add accessory letter to end of model number, i.e.: 3790C