

Lab-Line
Refrigerators/Freezers



Large Capacity Refrigerators



Model 3784

Our Largest Cold Storage Units!

- Storage of biological materials and other nonflammable substances
- Six models available—from 23, 49, and 72 cu. ft.
- CFC-free insulation and coolant
- ★ 60 Month Warranty on Compressor



Product Description

- Available in glass or solid door styles.
- Key locks on all doors.
- Adjustable thermostatic control from 1° to 10°C.
- Temperature uniformity of $\pm 0.5^{\circ}\text{C}$ at 4°C.
- .040 polished aluminum exterior.
- 22 gauge stainless steel floor.
- Foamed-in-place insulation.
- PVC coated adjustable shelves.
- Compressor bottom-mounted for easy access.
- UL listed.

PRODUCT SPECIFICATIONS

Model #	Interior Dimensions Inches (cm)			Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D	
3780 & 3780-1	23.75 (60)	61.5 (157)	26 (66)	27 (69)	78.25 (198)	31.75 (81)	245 (111)
3781 & 3781-1	23.75 (60)	61.5 (157)	26 (66)	27 (69)	78.25 (198)	31.75 (81)	277 (126)
3783 & 3783-1	50.5 (129)	61.5 (157)	26 (66)	54 (138)	78.25 (198)	31.75 (81)	415 (188)
3784 & 3784-1	50.5 (129)	61.5 (157)	26 (66)	54 (138)	78.25 (198)	31.75 (81)	480 (218)
3786 & 3786-1	74.75 (190)	61.5 (157)	26 (66)	78 (198)	78.25 (198)	31.75 (81)	585 (265)
3787 & 3787-1	74.75 (190)	61.5 (157)	26 (66)	78 (198)	78.25 (198)	31.75 (81)	655 (297)

*Exterior height does not include castors. Add 6" for castor height.

ORDERING INFORMATION

Model #	Volts	Electrical Hz	Capacity (cu. ft.)	Number of Shelves (sq. ft.)	Swinging Doors	Compressor	US List Price
3780 ¹	120	60	23	3 (13)	1 (Solid)	1/3 hp	\$4,047.00
3780-1CE ²	240	50	23	3 (13)	1 (Solid)	1/3 hp	\$4,047.00
3781 ¹	120	60	23	3 (13)	1 (Glass)	1/3 hp	\$4,598.00
3781-1CE ²	240	50	23	3 (13)	1 (Glass)	1/3 hp	\$4,598.00
3783 ¹	120	60	49	6 (29)	2 (Solid)	1/2 hp	\$5,398.00
3783-1CE ²	240	50	49	6 (29)	2 (Solid)	1/2 hp	\$5,398.00
3784 ¹	120	60	49	6 (29)	2 (Glass)	1/2 hp	\$6,342.00
3784-1CE ²	240	50	49	6 (29)	2 (Glass)	1/2 hp	\$6,342.00
3786 ¹	120	60	72	9 (43)	3 (Solid)	1/2 hp	\$7,260.00
3786-1CE ²	240	50	72	9 (43)	3 (Solid)	1/2 hp	\$7,260.00
3787 ¹	120	60	72	9 (43)	3 (Glass)	1/2 hp	\$8,677.00

¹ UL, ² CE

ACCESSORIES

Model #	Description	Model #	Description	US List Price
B*	Digital Temperature Display and Alarm with Remote Contacts	3780-8	Shelf for 3780, -1CE, and 3781, -1CE	\$142.00
L	6" Adjustable Stainless Steel Legs	3783-8	Shelf for 3783, -1CE, and 3784, -1CE	\$142.00
DT*	Digital Temperature Readout only	3786-8	Shelf for 3786, -1CE, and 3787, -1CE	\$142.00
S	Stainless Steel Interior and Exterior			

When ordering, add accessory letter to end of model number, i.e.:3782DT
Contact Customer Service for ordering information
* Not available on -CE models