

Lab-Line  
Refrigerators/Freezers



## Large Capacity Freezers



Model 3785



Model 3788

### Reach in for Easy Access!

- Storage of biological materials and other nonflammable substances
- Three models available—from 23, 49 and 72 cu. ft.
- CFC-free insulation and coolant
- ☆ 60 Month Warranty on Compressor



### Product Description

- Automatic defrost.
- Key locks on all doors.
- Adjustable thermostat control from -17°C to -23°C.
- Temperature control of ±0.5°C at -23°C.
- .040 polished aluminum exterior.
- 22 gauge stainless steel floor.
- Covered interior corners for ease of cleaning.
- Foamed-in-place insulation.
- PVC coated adjustable shelves.
- Compressor bottom-mounted for easy access.
- Swinging door(s).

### PRODUCT SPECIFICATIONS

Model #	Interior Dimensions Inches (cm)			Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D	
3782 & 3782-1CE	23.75 (60)	61.5 (156)	26 (66)	27 (69)	78.25 (199)	29.5 (75)	286 (130)
3785 & 3785-1CE	50.5 (128)	61.5 (156)	26 (66)	54 (137)	78.25 (199)	29.5 (75)	425 (193)
3788 & 3788-1	74.75 (190)	61.5 (156)	26 (66)	78 (198)	79.25 (201)	29.5 (75)	598 (271)

### ORDERING INFORMATION

Model #	Electrical Volts	Hz	Number of Capacity (cu. ft.)	Shelves (sq. ft.)	Solid Doors	Compressor	US List Price	Model #	Electrical Volts	Hz	Number of Capacity (cu. ft.)	Shelves (sq. ft.)	Solid Doors	Compressor	US List Price
3782 <sup>1</sup>	120	60	23	3 (13)	1	1/2 hp	\$5,100.00	3785-1CE <sup>2</sup>	240	50	49	6 (29)	2	3/4 hp	\$6,861.00
3782-1CE <sup>2</sup>	240	50	23	3 (13)	1	1/2 hp	\$5,100.00	3788 <sup>1</sup>	240	60	72	9 (43)	3	3/4 hp	\$9,512.00
3785 <sup>1</sup>	120	60	49	6 (29)	2	3/4 hp	\$6,861.00	3788-1	240	50	72	9 (43)	3	3/4 hp	\$9,512.00

<sup>1</sup> UL, <sup>2</sup> CE

### ACCESSORIES

Model #	Description	US List Price
3782-8	Shelf for 3782, -1CE, and 3782, -1CE	\$142.00
3785-8	Shelf for 3785, -1CE, and 3785, -1CE	\$140.00
3788-8	Shelf for 3788, -1CE, and 3788, -1CE	\$142.00