

APPLICATIONS



Pathology



General Laboratory



Laboratory Glassware

PRODUCTS



REVERSE OSMOSIS

SELECTION GUIDE

REVERSE OSMOSIS SYSTEMS

	Product	System Production	Product Water Flow ℓ /Hr 25°C/60 psi	Reservoir Capacity Options	Mounting	Page
The Ultimate RO System!	Diamond™ RO	From 50 to 400 ℓ /day	6-24	30, 60 or 100 liters	Wall or Bench	380
Smaller Size, But Not Smaller Volume!	EASYpure® RO	Up to 100 ℓ /day	10	30 liters	Wall or Bench	383
Larger Volume Requirements? No Problem!	ULTROpure	Up to 1325 ℓ /day	60, 80, 110, or 157	100+ liters	Wall or Bench	384
Simple and Efficient!	ROpure LP	From 100 to 400 ℓ /day	15, 20, 30, or 40	100+ liters	Wall	385
Type II Water Directly from Your Tap!	ROpure ST	From 50 to 350 ℓ /day	11	30 liters	Wall or Bench	386

STORAGE RESERVOIRS

	Product	Storage Capacity	Material of Construction	Mounting	Page
Compact Storage System!	Diamond RO Storage	60 liters	Fluorinated Polyethylene	Wall or Bench	382
Storage for All Your Water Requirements!	Diamond RO Storage	30 liters	Fluorinated Polyethylene	Wall or Bench	388
Our Most Versatile Reservoir!	100 Liter	100 liters	Polyethylene	Wall	389
Our Largest Storage Reservoirs!	Fiberglass Tanks	50, 100, Or 200 gallons	Premium Fiberglass Resin	Floor	390