



Diamond™ RO Storage Reservoirs



60 Liter Storage Reservoir



30 Liter Storage Reservoir

Compact Storage System!

- 30 and 60 liter storage capacity
- Automatic operation with Diamond RO
- Range of accessories

Product Description

- Complete with automatic controls, ensuring automatic operation when used in conjunction with ROpure Infinity® and DIamond RO reverse osmosis systems.
- Draw-off points include spigot, dishwasher port, and feed to NANOpure® DIamond systems.
- Constructed of rigid fluorinated polyethylene, providing strength and protection against material contamination of the water.
- Can be bench or wall mounted.
- Vent protected to eliminate the introduction of airborne contamination.
- Pump protector provided; protects NANOpure to ensure it will not be damaged when there is not enough water in the reservoir.
- 30 liter tank sends signal to DIamond RO displaying the water level on the RO control panel.

Optional Accessories

Sight Glass for 60 Liter Reservoir

- Available to alert you of water level within the reservoir.

Distribution Pump

- Perfect for supplying water to a remote NANOpure or glassware washer.
- Easily mounts to base of 60 liter Infinity reservoir.
- Can supply reverse osmosis water up to 35 feet away.
- Capable of delivering up to 1 gallon of water/min. at 15 psi.
- Automatically circulates water for 15 minutes every hour when in standby.
- Pump Protector: ensures distribution pump will not operate when tank is empty.

PRODUCT SPECIFICATIONS

Model #	Dimensions Inches (cm)			Shipping Weight Lb. (kg)	Operating Weight Lb. (kg)
	W	H	D		
D9021	20 (51)	24 (61)	17 (43.2)	35 (15.8)	168 (76)
D12681	20 (51)	18.3 (46.5)	15 (38)	23 (10.5)	104 (47)

ORDERING INFORMATION

Model #	Description	US List Price
D9021	60 Liter Storage Reservoir	\$978.00
D12681	30 Liter Storage Reservoir	\$750.00

ACCESSORIES

Model #	Description	US List Price	Model #	Description	US List Price
PU902X1	Optional Distribution Pump for D9021 60 Liter Reservoir, 120V	\$697.00	TU902X1	Optional Sight Glass for D9021 60 Liter Reservoir	\$92.60
PU902X2	Optional Distribution Pump for D9021 60 Liter Reservoir, 240V	\$697.00	AY1268X2	Optional Bench Stand for D12681 Reservoir	\$118.00