

Barnstead Reverse Osmosis



Storage Reservoirs



Diamond™ RO - 30 Liter Storage Reservoir - NANOpure® Diamond

Storage for All Your Water Requirements!

- 30 liters to 200 gallon capacity
- Strong, inert materials
- Completely automatic operation

Product Description

- Reverse osmosis systems require storage to eliminate back pressure on the reverse osmosis membrane and provide the quantity of water that will meet your requirements.
- Barnstead provides a variety of storage reservoirs ranging in volume from 30 liters to 200 gallons, depending upon your needs. All tanks are designed to ensure automatic operation when used in conjunction with Barnstead reverse osmosis systems.

CHOOSE THE REVERSE OSMOSIS SYSTEM AND STORAGE RESERVOIR TO MEET YOUR LABORATORY REQUIREMENTS

	30 Liter	60 Liter	100 Liter	Fiberglass Tanks 50, 100 & 200 gallon
EASypure® RO	Optional			
Diamond RO	Yes	Yes	Yes	Yes
ROpure™ LP			Yes	Yes
ULTROpure™			Yes	Yes

ORDERING INFORMATION

Model #	W	Dimensions Inches (cm) H	D	Operating Weight (Full) Lb. (kg)	Shipping Weight Lb. (kg)	US List Price
EASypure RO 30 Liter Storage Reservoir						
TY742X2A	18 (47.2)	29.5 (15)	14.7 (74.9)	130 (39)	61.5 (28)	\$1,297.00
ROpure Infinity and Diamond RO 60 Liter Storage Reservoir						
D9021	20 (51)	24 (61)	17 (43.2)	168 (76)	35 (15.8)	\$978.00
Diamond RO 30 Liter Storage Reservoir						
D12681	20 (51)	18.3 (46.5)	15 (38)	104 (47.2)	23 (10.5)	\$750.00

ACCESSORIES

Model #	Description	US List Price	Model #	Description	US List Price
D7424	Optional Bench Mounting Stand (to be used if bench mounting an EASypure RO when used in conjunction with an optional 30 Liter Tank)	\$174.00	TU902X1	Optional Sight Glass for 60 Liter Infinity Reservoir	\$92.60
PU902X1	Optional Distribution Pump for 60 Liter Reservoir, 120V	\$697.00	AY1268X2	Bench Mounting Stand for Diamond RO 30 Liter Reservoir	\$118.00
PU902X2	Optional Distribution Pump for 60 Liter Reservoir, 240V	\$697.00	16975	Float Assembly 50-200 Gallon Reservoirs	\$566.00