


Lab-Line Shakers



Dedicated Platforms with Clamps

ORDERING INFORMATION




Model 3524-3



Test Tube Platform



Model 3599-12



Model 3524-39



Model 3599-13



Model 3524-49

Lab-Line Shakers



ORDERING INFORMATION

18" x 18" (45.7 x 45.7 cm) Platforms with Clamps for
Shaker Series 3506, 4629, 3508, MaxQ 2000, MaxQ 4000

Model #	Vessels	Capacity	US List Price
3522-3	Adapter Set for Tiering 2 Platforms	—	\$61.00
3524-1*	10 ml Erlenmeyer Flask	113	\$682.00
3524-2*	25 ml Erlenmeyer Flask	64	\$586.00
3524-3*	50 ml Erlenmeyer Flask	64	\$586.00
3524-4*	125 ml Erlenmeyer Flask	32	\$586.00
3524-5*	250 ml Erlenmeyer Flask	25	\$360.00
3524-6*	300 ml Erlenmeyer Flask	16	\$360.00
3524-7*	500 ml Erlenmeyer Flask	16	\$417.00
3524-8	1 l Erlenmeyer Flask	9	\$417.00
3524-9	2 l Erlenmeyer Flask	5	\$367.00
3524-35*	10 ml Beaker	152	\$353.00
3524-36*	20 ml Beaker	104	\$353.00
3524-37*	30 ml Beaker	104	\$353.00
3524-38*	100 ml Beaker	44	\$353.00
3524-39*	250 ml Beaker	25	\$331.00
3524-40*	400 ml Beaker	25	\$331.00
3524-41*	600 ml Beaker	16	\$331.00
3524-42*	10 to 13 mm Tube Rack (4 Racks)	240	\$404.00
3524-48*	14 to 16 mm Tube Rack (4 Racks)	192 Tubes	\$404.00
3524-43*	18 to 20 mm Tube Rack (4 Racks)	176 Tubes	\$404.00
3524-47	16 to 25 mm Slant Tube Rack (2 Racks)	24 Tubes	\$405.00
3524-44*	22 to 25 mm Tube Rack (4 Racks)	96 Tubes	\$404.00
3524-45*	50 ml Centrifuge Tube Rack (4 Racks)	96 Tubes	\$404.00
3524-46*	Utility Tray	—	\$280.00
3524-49**	Separatory Funnel (500 ml to 2 l) (models 3506, 3508, MaxQ2000)	3	\$509.00

* Stackable

**For open air models only. Not for incubated shakers.

Dedicated Platforms: If using a single vessel size only, such as a flask, these platforms provide maximum capacity and come with clamps installed.

30" x 18" (76 x 45.7) Platforms with Clamps for
Shaker Series 3597, MaxQ 3000, MaxQ 5000

Model #	Vessel	Capacity	US List Price
3599-1	10 ml Erlenmeyer Flask	112	\$1,390.00
3599-2	25 ml Erlenmeyer Flask	112	\$1,390.00
3599-3	50 ml Erlenmeyer Flask	112	\$1,142.00
3599-4	125 ml Erlenmeyer Flask	56	\$1,094.00
3599-5	250 ml Erlenmeyer Flask	40	\$916.00
3599-6	300 ml Erlenmeyer Flask	24	\$899.00
3599-7	500 ml Erlenmeyer Flask	24	\$899.00
3599-8	1 L Erlenmeyer Flask	12	\$872.00
3599-9	2 L Erlenmeyer Flask	12	\$910.00
3599-10*	4 L Erlenmeyer Flask	6	\$865.00
3599-11*	6 L Erlenmeyer Flask	4	\$829.00
3599-12	2500 ml Low Form Culture Flask	5	\$1,085.00
3599-13	2800 ml Fernbach Flask	6	\$829.00
3599-42	10 to 13 mm Tube Rack (7 Racks)	420 Tubes	\$819.00
3599-48	14 to 16 mm Tube Rack (7 Racks)	336 Tubes	\$773.00
3599-43	18 to 20 mm Tube Rack (7 Racks)	308 Tubes	\$773.00
3599-44	22 to 25 mm Tube Rack (7 Racks)	168 Tubes	\$773.00
3599-45	50 ml Centrifuge Tube Rack (7 Racks)	168 Tubes	\$773.00
3599-46	Utility Tray	—	\$612.00
3599-49	Separatory Funnel (500ml to 2L)	10	\$1,513.00

36" x 24" (91 x 61 cm) Platforms with Clamps for Shaker Series MaxQ 3000

Model #	Vessel	Capacity	US List Price
3599-52	50 ml Erlenmeyer Flask	160	\$1,648.00
3599-53	125 ml Erlenmeyer Flask	80	\$1,648.00
3599-54	250 ml Erlenmeyer Flask	70	\$1,376.00
3599-55	300 ml Erlenmeyer Flask	70	\$1,262.00
3599-56	500 ml Erlenmeyer Flask	40	\$1,262.00
3599-57	1 L Erlenmeyer Flask	18	\$1,205.00
3599-58	2 L Erlenmeyer Flask	15	\$1,205.00
3599-59	4 L Erlenmeyer Flask	8	\$1,205.00
3599-60	6 L Erlenmeyer Flask	6	\$1,432.00
3599-61	2500 ml Low Form Culture Flask	6	\$1,432.00
3599-62	2800 ml Fernbach Flask	8	\$1,432.00
3599-42	10 to 13 mm Tube Rack (10 Racks)	600 Tubes	\$819.00
3599-48	14 to 16 mm Tube Rack (10 Racks)	480 Tubes	\$773.00
3599-43	18 to 20 mm Tube Rack (10 Racks)	440 Tubes	\$773.00
3599-44	22 to 25 mm Tube Rack (10 Racks)	240 Tubes	\$773.00
3599-45	50 ml Centrifuge Tube Rack (10 Racks)	240 Tubes	\$773.00
3599-96	Separatory Funnel (500 ml to 2 L)	16	\$2,472.00

* Not for use with 3525