

Lab-Line Shakers



Low Cost Orbital Shakers



Model 1345



Model

Economical Reliable Shakers!

- Compact design conserves bench space
- Adjustable speed control from 40 to 220 rpm
- Built-in timer for time dependent studies
- ☆ 24 Month Warranty

Product Description

- Solid state speed controller adjusts from 40 to 220 rpm.
- Operates in temperatures from 0°C to 40°C.
- Built-in timer can be set for continuous operation or for timed shaking up to 60 minutes.
- Durable powder coated cold rolled steel body construction.
- Model 1345 holds thirteen (13) 125 ml flasks; model 1346 holds nine (9) 250 ml flasks.
- Shakers include platform and clamps.

PRODUCT SPECIFICATIONS

Model #	Flask Size (ml)	Platform Capacity	Platform Dimensions Inches (cm) W x D		Overall Dimensions Inches (cm) W x H x D			Shipping Weight Lb (kg)
			W	D	W	H	D	
1345 ¹	125	13	12.25 (31)	12.25 (31)	12.25 (31)	6.25 (16)	13.4 (34)	18 (8)
1345-1	125	13	12.25 (31)	12.25 (31)	12.25 (31)	6.25 (16)	13.4 (34)	18 (8)
1346 ¹	250	9	12.25 (31)	12.25 (31)	12.25 (31)	6.25 (16)	13.4 (34)	18 (8)
1346-1	250	9	12.25 (31)	12.25 (31)	12.25 (31)	6.25 (16)	13.4 (34)	18 (8)

¹ ETL

ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)		Maximum Load Capacity Lbs (kg)	Orbit Size Inches (cm)	Speed Range (rpm)	US List Price
		Amps	Watts				
1345 ¹	120	0.4	50	10 (5)	0.75 (1.9)	40-220	\$238.40
1345-1	240	0.2	50	10 (5)	0.75 (1.9)	40-220	\$844.00
1346 ¹	120	0.4	50	10 (5)	0.75 (1.9)	40-220	\$844.00
1346-1	240	0.2	50	10 (5)	0.75 (1.9)	40-220	\$844.00

Operating Environment: 0-40°C, 20-80% RH, non-condensing