



Reciprocating Shaker



Model 3506

Product Description

- Variable speed control from 10 to 250 rpm.
- Provides linear shaking with an adjustable stroke length of 0.5" to 1.5" (1.2 to 3.8 cm)
- Interchangeable 13" x 11" (33 x 28 cm) or 18" x 18" (46 x 46 cm) or five tiered 11" x 13" (28 x 33 cm) shaking platforms are available

Versatile Reciprocating Operation!

- Variable speed control from 10 to 250 rpm
- Adjustable shaking motion
- Built-in timer for repeatable results
- ☆ 24 Month Warranty

to accommodate a variety of laboratory vessels and glassware. (See Accessories)

- Built-in timer can be set from 1 minute to 60 minutes or for continuous operation.
- Durable powder coated steel outer body.
- Operates in temperatures from 4°C to 40°C and 20 to 80% humidity conditions non-condensing.
- Load capacity 15 lb. (6.8 kg).

ORDERING INFORMATION

Model #	Volts	Amps	Watts	W	Overall Dimensions Inches (cm)			Shipping Weight Lbs. (kg)	US List Price
					H	D			
3506	120	1.0	120	12.2 (31)	7.4 (19)	17.1 (43)	37(17)	\$96.50	
3506-1	240	0.5	120	12.2 (31)	7.4 (19)	17.1 (43)	37(17)	\$1,634.00	

Operating Environment: 4-40°C; 20-80% RH, non-condensing Speed Range: 10-250 rpm

ACCESSORIES

Model #	Description	Platform Size Inches (cm)	US List Price	Model #	Description	Platform Size Inches (cm)	US List Price
30100	Universal Platform	13"x11" (33x28)	\$141.00	3523-39	Dedicated separatory funnel platform 500ml to 2 L (holds 2 funnels)	13" x 11" (33x28)	\$319.00
30110	Universal Platform	18"x18" (46x46)	\$214.00	3524-49	Dedicated separatory funnel platform 500ml to 2L (holds 3 funnels)	18" x 18" (46x46)	\$509.00
30152	Universal flask clamp - 50 ml		\$18.00	3523-51	Five Level, Tiered Blood Bag Platform	13" x 11" x 13.8" (33x28x35)	\$280.00
30153	Universal flask clamp - 125 ml		\$18.00				
30154	Universal flask clamp - 250 ml		\$21.00				
30156	Universal flask clamp - 500 ml		\$24.00				
30157	Universal flask clamp - 1000 ml		\$33.00				



Model 3523-1

Product Description

- A five tier aluminum platform assembly measuring 11" x 13" (28 x 33 cm).
- Distance between platforms is 3.25" (8.26 cm.)
- Each of the five platforms is equipped with a white ribbed vinyl mat.
- Mixing contents of 10" x 10" (25.4 x 25.4 cm) bags of platelets used in Platelet Pheresis applications.

Applications

- Platelet pheresis
- Gel staining
- Separatory funnel protocols

ORDERING INFORMATION

Model #	W	Overall Dimensions Inches (cm)		D	Shipping Weight Lb. (kg)	US List Price
		H				
3523-51	11 (28)	13.8 (34.9)		13 (33)	11 (5)	\$280.00