

Lab-Line Shakers



Dual Action Orbital/Reciprocating Shaker



Model 3508

**For a complete system order:
1. Shaker 2. Platform 3. Clamps**

Product Description

- Provides 0.4" (0.9 cm) orbital or linear shaking.
- Variable speed control accurately regulates shaker speed from 50 to 400 rpm controlled by a rotary dial for both shaking motions.
- Control lever changes shaking motion while unit is stationary or in operation.
- Load capacity 15 lb. (7 kg).
- A variety of 13" x 11" (38 x 28 cm) or 18" x 18" (46 x 46 cm).

Two Mixing Actions— One Low Price!

- Provides both reciprocating and orbital motion
- Variable speed from 50 to 400 rpm
- Shaker operates in continuous timer mode
- Molded polyethylene body - won't chip, rust or corrode
- ☆ 24 Month Warranty

interchangeable platforms are available to accommodate many laboratory vessels and glassware.

- Units operate in temperature ranges of 0 to 40°C and humidity conditions from 20 to 80% non-condensing enabling them to be used in incubators, warm rooms and environmental chambers.

Applications

- Bacteria culture
- Staining and destaining gels
- Blotting
- Platelet pheresis applications

ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)		W	Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)	US List Price
		Amps	Watts		H	D		
3508	120	0.8	100	11.5 (29)	7.75 (19.7)	14 (36)	36 (16)	\$1,623.00
3508-1	240	0.4	100	11.5 (29)	7.75 (19.7)	14 (36)	36 (16)	\$1,623.00

Speed range: 50 to 400 rpm

Orbit Diameter: 0.4" (0.9cm)

Max Load: 15Lbs (7kg)

Operating Environment: 0 to 40°C, 20 to 80% RH, non-condensing

ACCESSORIES

Universal Platforms			Clamps for Universal Platforms				
Model #	Description	US List Price	Model #	Vessel Size	Maximum Number of vessels per platform size (Inches)		US List Price
					11x13	18x18	
30100	Universal 11"x13" (27.9x33cm)	\$141.00	30175	Microplate/Deep Well Plate	5	10	\$18.00
30110	Universal 18"x18" (45.7x45.7cm)	\$214.00	30150	10 ml Flask	60	113	\$16.00
			30151	25 ml Flask	60	64	\$16.00
			30152	50 ml Flask	30	32	\$18.00
			30153	125 ml Flask	13	32	\$18.00
			30154	250 ml Flask	9	16	\$21.00
			30156	500 ml Flask	6	16	\$24.00
			30157	1 L Flask	4	5	\$33.00